

## PATIENT IMMOBILIZER CLEANING INSTRUCTIONS

The TZ Medical **Patient Immobilizer, GE Patient Immobilizer, and XXL Extender** are made of strong, durable vinyl material that will withstand many uses; however, care must be taken when cleaning the products to keep them in good condition.

**Please note:** *These cleaning instructions are only for the Patient Immobilizer, GE Patient Immobilizer, and XXL Extender; the wrist wraps are indicated for **single use only**.*

### Washing:

1. **Prewash** to remove any caked-on material by rinsing with water and wiping with a non-abrasive cloth or sponge.
2. **Wash** the Patient Immobilizer with a mild detergent (such as ENZOL<sup>®</sup> Fast-acting Enzymatic Cleaner) by briefly spraying and wiping down. Isopropyl alcohol (rubbing alcohol) may be used to spot clean only but **NOT** for soaking. Avoid abrasives such as brushes or scouring pads as they may scratch the surface.
3. **Rinse** off all cleaning agents. *Only clean the Patient Immobilizer with chemicals which are indicated for cleaning vinyl surfaces.*
4. **Dry** using a cloth, air, or heat up to 120 °F.

### CAUTION:

- **DO NOT allow the Patient Immobilizer to be heated above 120 °F.** If a Patient Immobilizer is suspected to have been heated above 120 °F, inspect for visible damage (melted, deformed, or discolored areas).
- **DO NOT steam clean** or allow the Patient Immobilizer to come into contact with steam as it can degrade the materials of the product.
- **DO NOT clean the Patient Immobilizer with any chemical not specifically indicated for the cleaning of vinyl surfaces**
- **DO NOT use a Patient Immobilizer that has evidence of damage to any part of the product.**