



TZ Medical
Sparked by your ideas

Panni XL
PANNICULUS RETRACTOR



1 - 4" x 14" Pad
1 - 12.5" x 13" Strap
(with lower adhesive strip)

The Panni® XL Retractor is a sterile pannus management system designed to retract and reposition excess tissue away from the incision site during and after C-section procedures. The Panni XL can be quickly applied by one person and is completely adjustable to ensure optimal patient comfort. It is shower compatible and can be worn for up to 3 days postpartum to assist with wound health.



Part#: CZ-PPR-XL

Features

- Designed for surgical and postpartum use
- Sterile
- Easily applied by one person, in seconds
- Completely adjustable
- Residue and latex-free
- Can be worn while showering
- Can be worn for up to 3 days



FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

21814 R0.2