



**TZ Medical**  
Sparked by your ideas


# PUSH-CLICK Y Connector Kit

The TZ Medical Push-Click Y Connector is designed to provide rapid single-handed manipulation of the valve with a simple push-and-click action. When the mechanism is compressed, the release lever holds the valve in an “open” position, allowing for free movement of wires and catheters. Press the release lever to simply and rapidly close the valve. The Push-Click Y Connector kit includes a Wire Introducer and Torque device.



## Features

- Simple one-handed operation
- Rapid and robust valve closure with audible actuation



Item Number	Description	UOM
-------------	-------------	-----

IDD-0030	Push-Click Y Connector kit with Wire Introducer and Torque Device	Box of 5
----------	---	----------

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

22256 R0.2