



The TZ Medical URPro is a universal system that directly replaces your current MAVIG™ shields with high-quality, long-lasting radiation protection panels. The new URPro Adapter panels are uniquely designed to fit over any MAVIG™ frame system* using an adjustable strap that quickly secures around the frame, making placement or removal quick and simple. More importantly, each panel is well-structured, with fabric binding around the edges to reduce common issues with fraying and deterioration caused by long-term usage.

**TZ Medical only supplies URPro Radiation Protection Shields/Panels. All other components are manufactured and trademarked products of the MAVIG™ corporation.*

Features

- Attenuates 98% of secondary and scatter radiation at 80 kVp
- Lead equivalent of 0.50 mm
- Flexible straps that help secure the panels to current MAVIG™ bedside drape framing systems
- Can replace any number of damaged panels
- Fabric binding and structured frame reduce long-term fraying and deterioration
- Available in 30" and 35" length options



Item Number	Description	UOM
BDR-1015	URPro Adapter Single 35" Panel	1 Each
BDR-1015-5	URPro Adapter 35" Panel Set	Set of 5
BDR-1025	URPro Adapter Single 30" Panel	1 Each
BDR-1025-5	URPro Adapter 30" Panel Set	Set of 5

FIND OUT MORE

