




The TZ Medical URPro is a universal system that directly replaces your current MAVIG™ shields with high-quality and long-lasting radiation protection panels. The URPro panels are uniquely designed to match your current MAVIG™ frame system*, allowing you to replace any and/or all damaged panels through a simple installation process. More importantly, each panel is well-structured with fabric binding around the edges to reduce common issues with fraying and deterioration caused by long-term usage.

*TZ Medical only supplies URPro Radiation Protection Shields/Panels. All other components are manufactured and trademarked products of the MAVIG™ corporation.

Features

- Attenuates 98% of secondary and scatter radiation at 80 kVp
- Lead equivalent of 0.50 mm
- Enhanced lower body radiation protection
- Directly compatible with current MAVIG™ bedside drape framing systems
- Simple installation process
- Can replace any number of damaged panels
- Fabric binding and structured frame reduces long-term fraying and deterioration
- Available in 30" and 35" length options



OLD PANELS

NEW PANELS

Item Number	Description	UOM
BDR-1002	URPro 35" Replacement Panel	1 panel
BDR-1002-5	URPro 35" Replacement Set	5 panels
BDR-1020	URPro 30" Replacement Panel	1 panel
BDR-1020-5	URPro 30" Replacement Set	5 panels

FIND OUT MORE

