




TZ Medical
Sparked by your ideas

SR HEADBAND



*The TZ Medical SR Headband
has a 0.25 mm lead equivalency*

When it comes to radiation, protecting your brain is a top priority. The TZ Medical SR Headband is a cost-effective way to block scatter radiation from reaching one of the body's most sensitive areas during diagnostic and interventional procedures.

The SR Headband is disposable, lead-free, easily worn, adjustable to any head size, and attenuates 87% of radiation at 80 kVp. It's a must-have for physicians, techs, and anyone involved in imaging procedures.

Features

- Attenuates 87% of scatter radiation at 80 kVp
- Comfortable and lightweight
- Lead-free and disposable
- Fully adjustable and can be worn over any style of cap
- Cost-effective way to increase safety in the lab



Item Number	Description	UOM
CZ-013	SR Headband - radiation protection headband	Box of 10

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

14115 R1.2