

# ACCUFLATE INFLATION DEVICE


The TZ Medical AccuFlate is a 20ml disposable device designed to accurately control and monitor pressure during inflation and deflation of dilation balloon catheters. Our patented technology allows the Inflation Device to generate pressure up to 30 atm or 40 atm with accuracy of  $\pm 1$  atm.



**Available in 30 atm  
and 40 atm options**



**Patented self-locking mechanism for easy one-handed locking control, as well as gradual and rapid pressure release ergonomic handle.**



**Precision manometer with 30 atm and 40 atm maximum pressure options.**



**Optional, male rotating Luer Lock, 3-way stopcock with "ON" or "OFF" tap handle.**

Item Number	Description	UOM
IDD-0001	30 atm Balloon Inflation Device	Box of 5
IDD-0040	40 atm Balloon Inflation Device	Box of 5
IDD-0010	30 atm AccuFlate with PTCA Kit	Box of 5 sets
IDD-0020	PTCA Accessory Kit	1 each



**PTCA Accessory Kit also available**

**FIND OUT MORE**

Patent: CN103028153A



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

19337 R1.4