



TZ Medical
Sparked by your ideas



OPERATING ROOM

PRODUCT SOLUTIONS

Improving healthcare through innovative product solutions that meet and exceed the expectations of our customers.



tzmedical.com



800.944.0187



info@tzmedical.com



TZ Medical
Sparked by your ideas

TABLE OF CONTENTS

| | | |
|------------------------------|-------|----|
| Patient Positioning Products | _____ | 4 |
| Lab/Hospital Essentials | _____ | 10 |
| Electrodes & Grounding Pads | _____ | 14 |
| Neurosurgery Products | _____ | 16 |

TZ Medical Inc.

20497 SW Teton Ave, Suite A

Tualatin, OR 97062

Office: 800.944.0187 Fax: 503.639.0239

www.tzmedical.com

POSITIONING PRODUCTS

PANNUS RETENTION SYSTEM



The TZ Medical Pannus Retention System is a hook and loop pad and strap system designed to safely and effectively retain excess tissue for high BMI patients. It can be used as a patient positioning tool for various procedural tasks, including gaining vascular groin access and unilateral hip replacements.

- **Hypoallergenic, latex and residue-free**
- **Adhesive pads are designed to reduce damage to the skin**
- **Ideal for groin access, anterior hip procedures, and more**
- **Completely adjustable system - will fit any sized patient**
- **It can be applied in less than one minute**
- **Easy to reposition after application**

Safely retains large amounts of tissue



Intended for skin applications, the proprietary adhesive on the pads is residue-free. It minimizes the risk of skin tears and degradation typically seen with tape and competing products.

| Item Number | Description | UOM |
|------------------|--|------------------|
| CZ-PRS-02 | Non-Sterile Two Pad and One Strap System | Box of 10 |
| CZ-PRS-04 | Sterile Four Pad and Two Strap System | Box of 10 |



CERVICAL & LUMBAR SPINE RETENTION SYSTEM



The PRS™ Cervical Spine Retention System is a simple and safe solution to position the patient for optimal access and better imaging by deflecting and positioning the patient's shoulders down towards the foot of the bed.

The PRS™ Lumbar Spine Retention System provides an optimal solution for safely positioning patients during lateral spine procedures by creating three anchor points on the patient and both sides of the procedure table, securing the patient comfortably on their side.

- Eliminates risk of joint injury caused by pulling wrist restraints for deflection
- Device will not release on its own during the procedure
- Adhesive pads are designed to reduce damage to the skin



| Item Number | Description | UOM |
|------------------|--|------------------|
| CZ-PRS-06 | Cervical Spine Retention System - Four pads, Two Straps, Non-Sterile | Box of 10 |
| CZ-PRS-09 | Lumbar Spine Retention System - Nine pads, Three Straps, Non-Sterile | Box of 10 |


POSITIONING PRODUCTS

CRANIAL IMMOBILIZATION SYSTEM



The Cranial Immobilization System (PRS-03) is a safe and effective solution for immobilizing the patient's head during imaging and various surgical procedures. This unique system, designed with a micro hook and loop technology, is completely adjustable to ensure optimal patient support and comfort.

- Immobilize the patient's head in a safe and effective manner
- Hypo-allergenic, latex-free and residue-free
- Completely adjustable system - easy to apply and reposition if needed
- Prolong equipment life-span and increase regulatory compliance by avoiding tape and tape residue



Great for Breast Deflection

Simple, effective breast deflection for pacemaker and ICD implantation; PRS-03 securely and safely deflects breast tissue away from the implant location, ensuring optimal accessibility.

| Item Number | Description | UOM |
|-------------|-------------------------------|-----------|
| CZ-PRS-03 | Cranial Immobilization System | Box of 10 |



STERILE DRAPE HOLDER



The TZ Medical Sterile Drape Holder is an adaptable solution that slides under any mattress and provides simple drape positioning for the comfort of both the patient and medical staff. Crafted with a durable polypropylene base and two adjustable 30" flexible goosenecks, the Drape Holder offers versatility in shaping to optimally position your patients and ensure their comfort.

- Eliminates the need for using bulky, large-footprint IV poles for drape positioning
- Reduces equipment congestion at the head of the bed and creates more room for medical staff
- Flexible gooseneck arms can be molded to any desired shape or height to accommodate patients of all shapes and sizes
- Compatible with all bed types; no arm rails or additional equipment necessary
- Sturdy, easy-to-clean polypropylene base
- Easily repositioned to adapt to changing procedural needs



| Item Number | Description | UOM |
|-------------|----------------------|------|
| SDH-0005 | Sterile Drape Holder | Each |

POSITIONING PRODUCTS

SOFTGEL POSITIONERS



TZ-GELDA



TZ-GELDA-3



TZ-GELAF



TZ-GELAC



TZ-GELPD



TZ-GELND



TZ-GELUPL



TZ-GELUPS



TZ-GELHRS



TZ-GELPHR



TZ-GELCC



TZ-GELHC



TZ-GELSHP



TZ-GELHDP



TZ-GELHP



TZ-GELAS



TZ-GELPP10



TZ-GELCR4-20



TZ-GELPP14



TZ-GELCR6-12



TZ-GELPP18



TZ-GELCR6-16



TZ-GELCR4-12



TZ-GELCR6-20



TZ-GELCR4-16



TZ-GELFT



TZ-GELFT3



TZ-GELFT3-FT



| Item Number | Description - Head Gels | UOM |
|-------------------|---|---------------|
| TZ-GELDA | Adult Gel Donut 8"OD x 3"ID x 1.75"H | 1 Each |
| TZ-GELDA-3 | Extra Deep Donut 8"OD x 3"ID x 3"H | 1 Each |
| TZ-GELPD | Pediatric Donut 5.5"OD x 2.25"ID x 1.25"H | 1 Each |
| TZ-GELND | Neonatal Donut 3.07"OD x 1.68"ID x 0.81"H | 1 Each |
| TZ-GELHRS | Adult Horseshoe Donut 8"OD x 3"ID x 1.75"H | 1 Each |
| TZ-GELPHR | Pediatric Horseshoe Donut 5.5"OD x 2.25"ID x 1.25"H | 1 Each |
| TZ-GELSHP | Supine Headrest 11"OD x 9.5"ID x 2.625"H | 1 Each |
| TZ-GELHDP | Head Pad 11"OD x 9.5"ID x 3.5"H | 1 Each |

| Item Number | Description - Body Gels | UOM |
|---------------------|--|---------------|
| TZ-GELPP10 | Gel Patient Positioner 10" x 4.5" x 2" (Sandbag Style) | 1 Each |
| TZ-GELPP14 | Gel Patient Positioner 14" x 4.5" x 2" (Sandbag Style) | 1 Each |
| TZ-GELPP18 | Gel Patient Positioner 18" x 4.5" x 2" (Sandbag Style) | 1 Each |
| TZ-GELCR4-12 | Flat Bottom Gel Chest Roll 12" x 4" x 4" | 1 Each |
| TZ-GELCR4-16 | Flat Bottom Gel Chest Roll 16" x 4" x 4" | 1 Each |
| TZ-GELCR4-20 | Flat Bottom Gel Chest Roll 20" x 4" x 4" | 1 Each |
| TZ-GELCR6-12 | Flat Bottom Gel Chest Roll 12" x 6" x 6" | 1 Each |
| TZ-GELCR6-16 | Flat Bottom Gel Chest Roll 16" x 6" x 6" | 1 Each |
| TZ-GELCR6-20 | Flat Bottom Gel Chest Roll 20" x 6" x 6" | 1 Each |

| Item Number | Description - Arm & Foot/Ankle Gels | UOM |
|------------------|--|---------------|
| TZ-GELAF | Flat Armboard Gel 24" x 4.5" x 0.625" | 1 Each |
| TZ-GELAC | Gel Contoured Armboard 20" x 5.75" x 2.25" | 1 Each |
| TZ-GELUPL | Ulnar Nerve Protector, Large 20.13" x 13.75" x 0.25" | 1 Each |
| TZ-GELUPS | Ulnar Nerve Protector, Small 16.25" x 6.5" x .5" | 1 Each |
| TZ-GELCC | Gel Calf Cup | 1 Each |
| TZ-GELHC | Gel Heel Cup | 1 Each |
| TZ-GELHP | Heel Protector Pad | 1 Each |
| TZ-GELAS | Ankle Stirrup Pad | 1 Each |

| Item Number | Description - Table Gels | UOM |
|---------------------|--|---------------|
| TZ-GELFT | Full Length Table Gel Overlay 20" x 72" x 0.5" | 1 Each |
| TZ-GELFT3-FT | Gel Foot Table Overlay 20" x 20" x 0.5" | 1 Each |
| TZ-GELFT3-HD | Gel Head Table Overlay 20" x 10" x 0.5" | 1 Each |
| TZ-GELFT3-T | Gel Torso Table Overlay 20" x 42" x 0.5" | 1 Each |
| TZ-GELFT3 | 3 Piece Table Gel Overlay (Includes: TZ-GELFT3-FT, TZGELFT3-HD, TZ-GELFT3-T) | 1 Each |

LAB/HOSPITAL ESSENTIALS

OPERATING ROOM MATTRESSES

Multi-layered, high-resiliency foam for added comfort, support, and durability. Designed for shorter procedures or short stay procedural/exam tables.

- Standard and memory foam mattress options available
- Sealed seams prevent liquid penetration
- Anti-bacterial, anti-fungal, and water/stain resistant
- Secondary fluid barrier protects foam if the cover is compromised
- Custom-sizing available for all mattresses and arm board pads
- Made in the USA



INDUSTRY-LEADING COMFORT

Our pressure-mapped custom foam layering ensures an optimal level of comfort for every procedure.

EASILY CLEANED REMOVABLE COVERS

Our expertly made covers provide a safe, comfortable, and durable solution for the healthcare industry.

DURABILITY YOU CAN TRUST


Advanced stability layers provide long-lasting comfort without sagging, ensuring consistent, hard-working performance every time it's used.



ANTI-FATIGUE MATS

OR Anti-Fatigue Mats are scientifically engineered with conforming gel and supportive foam to provide the ultimate in comfort and support to surgeons, nurses, and technicians. Mats are water-resistant, easy to clean, and provide the proper friction level for damp or dry surgical booties.

- Reduces and prevents physical discomfort and fatigue
- Properly conforms and supports feet while promoting proper circulation
- Relieves pressure and restores energy while standing
- Improve comfort and ergonomics during procedures
- Reduce the number of stretches and postural changes



| Item Number | Description | UOM |
|-----------------|--|---------------|
| GP-18X24 | OR Durable 3/4" Thick Bio-Foam® Size 18X24, Blue | 1 Each |
| GP-20X32 | OR Durable 3/4" Thick Bio-Foam® Size 20X32, Blue | 1 Each |
| GP-20X48 | OR Durable 3/4" Thick Bio-Foam® Size 20X48, Blue | 1 Each |
| GP-30X72 | OR Durable 3/4" Thick Bio-Foam® Size 30X72, Blue | 1 Each |

LAB/HOSPITAL ESSENTIALS

PRO11 HEADSET SYSTEM



The TZ Medical Pro11 Headset System enables teams to talk hands-free in the operating room and adjacent rooms - no raised voices or straining to hear, just clear, calm communication between every member of the team, reducing stress and creating a better working environment. The system has full connectivity with external broadcast and audiovisual streaming facilities, enabling the physician to participate real-time in remote broadcasts and seminars. Using high-quality DECT frequencies, multiple systems can be used in the same premises without interference.



| Item Number | Description | UOM |
|----------------------|---|-------------------|
| TZ-Q-P11BS | Pro11 Base Station | 1 Each |
| TZ-Q-P11CH | Pro11 Charger | 1 Each |
| TZ-Q-P11HS | Pro11 Headset (battery & headband included) | 1 Each |
| TZ-Q-P11BAT | Battery for Headset | 1 Each |
| TZ-Q-P11HB | Headband for Headset | 1 Each |
| TZ-Q-P11NB | Neckband for Headset | 1 Each |
| TZ-Q-RH7HC | Hygiene Covers | Box of 100 |
| TZ-Q-RH7EC | Ear Cushions | Box of 10 |
| TZ-Q-P10SHELF | Battery Charger Shelf | 1 Each |

STAINLESS STEEL INSTRUMENT TABLES

TZ Medical now offers a line of stainless steel instrument tables. Each lightweight stainless steel table has heavy-duty swivel casters, a lower shelf, fully welded legs, and 16-gauge steel. Available in three lengths, TZ tables provide a stable, easy-to-clean surface for all procedures throughout the hospital.



- 16-gauge stainless steel
- Closed hems for simple sterilization
- Convenient lower shelf for storage
- Comes in three lengths: 48", 60", and 72"
- Made in the USA



| Item Number | Description | UOM |
|----------------|--|---------------|
| SSIT-48 | TZ Medical Instrument Table - 48"L x 24"W x 34"H | 1 Each |
| SSIT-60 | TZ Medical Instrument Table - 60"L x 24"W x 34"H | 1 Each |
| SSIT-72 | TZ Medical Instrument Table - 72"L x 24"W x 34"H | 1 Each |

ELECTRODES & GROUNDING PADS

STERILE DEFIBRILLATION ELECTRODES



Our multi-function sterile electrodes are ideal for a variety of procedures, including TAVR, AVR, MiniMaze, ICD Implants, Laser Lead Extraction and much more. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Our sterile electrodes offer direct connection to Zoll® defibrillators. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.

- Lead wire lengths available in 60" and 120" for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Adult, pediatric and extra large backpad options available
- Radiotransparent and radiotranslucent options available

| Description | Size | LifePak® ¹ | Universal ² | Zoll® ³ |
|---|------------|---|------------------------|--|
| Adult Radiotranslucent, Sterile | A | P-121-T1 <i>PAD 1111 Adapter*</i> | P-121-T1 | P-121-T1 <i>PAD 101PD Adapter*</i> |
| Adult Radiotransparent Pads and Leadwires, Large Posterior Pad, Sterile | A/B | P-224-T1 <i>PAD 1111 Adapter*</i> | P-224-T1 | P-224-Z1 |
| Pediatric Radiotranslucent Pads, Sterile | C | P-122-T1 <i>PAD 1111 Adapter*</i> | P-122-T1 | P-122-T1 <i>PAD 101PD Adapter*</i> |

Pad Dimensions:

A = 6"x4.25"

B = 6"x9"

C = 4.25"x2.876"

**Adapters are reusable*

1 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)

2 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

Pediatric Use: Pediatric patients whose weight is less than 25Kg (55 lbs).



GROUNDING PADS

TZ Medical's Grounding Pads are indifferent skin contact electrodes designed for energy return of radiofrequency energy in surgical and electrophysiology procedures. The soft foam structure and sticky hydrogel provide optimal adhesion that comfortably conforms to the patient's body while maintaining a large surface area to focus radiofrequency energy at the application site.

- Utilized for safe return of electrosurgical energy from the patient
- 9.8 foot Lead Wire
- Designed with adhesive hydrogel for secure patient application
- Compatible with most generators on the market
- Comfortable foam structure conforms to the patient's body



| Item Number | Description | UOM |
|-------------|-------------------------------------|-----------|
| P-411-V1 | Disposable Grounding Pad Electrodes | Box of 50 |

NEUROSURGERY PRODUCTS

DISTRACTION PINS



Thread Comparison




TZ Medical
Titanium
Distraction Pin



Stainless Steel
Distraction Pin

TZ Medical titanium distraction pins for anterior cervical fusions have a self-drilling design that allows for effortless insertion. In addition, TZ Medical's distraction pins have a 10% wider diameter than typical distraction pins for increased durability and complete surgical confidence. TZ Medical distraction pins are color-anodized and laser-engraved for easy visual identification of size. They are available to order in four sizes: 12mm, 14mm, 16mm, and 18mm.

- Triple-fluted design for easier insertion into bone
- Wider diameter than typical distraction pins
- Color-coded sizes for faster identification
- Made with 100% titanium
- Cost-effective
- Safer and sturdier than pins made from stainless steel
- Compatible with all current ACF systems



| Item Number | Description | UOM |
|------------------|--------------------------------------|------------------|
| DP-12-TB | 12mm Titanium Distraction Pin Blue | Box of 10 |
| DP-14-TY | 14mm Titanium Distraction Pin Yellow | Box of 10 |
| DP-16-TG | 16mm Titanium Distraction Pin Green | Box of 10 |
| DP-18-TP | 18mm Titanium Distraction Pin Purple | Box of 10 |
| DP-S | Caspar Screwdriver Universal | 1 Each |
| DP-72-CDR | Caspar Distractor Right | 1 Each |
| DP-74-CDL | Caspar Distractor Left | 1 Each |



SKULL PINS



TZ Medical Skull Pins are a disposable solution for cranial stabilization. The stainless steel skull pins are inserted into the skull clamp to fixate the patient's head during surgery. The pins are packaged sterile in a pack of three and come with a protective cap, making them safe and easy to handle.

- Cost-effective
- Improves cranial stabilization
- Safe and easy to handle
- Sharp and durable stainless steel
- Works with existing skull clamps



| Item Number | Description | UOM |
|---------------|--|------------------|
| SP-001 | Adult, Stainless Steel Skull Pins (3 pins per package) | Box of 10 |



TZ Medical

Sparked by your ideas

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com