



**TZ Medical**  
Sparked by your ideas


# ADULT DEFIB PADS



Our multi-function electrodes are designed for EP and Cath Labs, Cardiac Surgery, Critical-Care Areas and other related departments. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Compatible with biphasic defibrillators, our electrodes offer direct connections to a variety of defibrillators, such as Zoll® and Medtronic® QUIK-COMBO™. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.



## Radiotransparency



## Features

- 60" lead wires for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Radiotransparent and radiotranslucent options available

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)




800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

10943 R3.8

# Pad Placement Guide



## Lead Wire Length

**TZ Medical: 60"**

**Competitor Electrodes: 42"**

**Competitor Electrodes: 24"**

Medtronic® M1 Connector



Universal T1 Connector



Zoll® Z1 Connector



Product Description	Size	Medtronic® <sup>1</sup>	Universal <sup>2</sup>	Zoll® <sup>3</sup>
Adult Radiotranslucent pads, leads out	A	P-111-M1	P-111-T1	P-111-Z1
Adult Radiotransparent pads and leadwires	A	P-211-M1	P-211-T1	P-211-Z1
Adult Radiotransparent pads and leadwires, leads out	A	P-311-M1		P-311-Z1
Adult Radiotransparent, large posterior pad	A/B	P-214-M1	P-214-T1	P-214-Z1

### Pad Dimensions:

**A=6"x 4.25"**

**B=6"x9"**

1 – M1 Medtronic models are comparable to Medtronic® Physio-Control QUIK-COMBO™ electrodes.

2 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)

3 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

10943 R3.8