



TZ Medical
Sparked by your ideas

ADULT DEFIB PADS




TZ Medical Defibrillation Electrodes are constructed to adhere and conform to the patient's body, providing superior current distribution, which results in better clinical outcomes.

Our multi-function electrodes are designed for EP and Cath Labs, Cardiac Surgery, Critical-Care Areas and other related departments. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Compatible with biphasic and monophasic defibrillators, our electrodes offer direct connections to a variety of defibrillators, such as Zoll® and Medtronic® QUIK-COMBO™. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.



Radiotransparency



Features

- 60" lead wires for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Radiotransparent and radiotranslucent options available

FIND OUT MORE



tzmedical.com




800.944.0187




info@tzmedical.com

10943 R3.6

Pad Placement Guide



Lead Wire Length



Medtronic® M1 Connector



Universal T1 Connector



Zoll® Z1 Connector



Product Description	Size	Medtronic® ¹	Universal ²	Zoll® ³
Adult Radiotranslucent pads, leads out	A	P-111-M1	P-111-T1	P-111-Z1
Adult Radiotransparent pads and leadwires	A	P-211-M1	P-211-T1	P-211-Z1
Adult Radiotransparent pads and leadwires, leads out	A	P-311-M1		P-311-Z1
Adult Radiotransparent, large posterior pad	A/B	P-214-M1	P-214-T1	P-214-Z1

Pad Dimensions:

A=6" x 4.25"

B=6" x 9"

1 – M1 Medtronic models are comparable to Medtronic® Physio-Control QUIK-COMBO™ electrodes.
 2 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)
 3 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

10943 R3.6