



TZ Medical
Sparked by your ideas

CERVICAL SPINE RETENTION SYSTEM



Four 5.5" x 7.5" Pads



Two 96" Straps

4 Pad, 2 Strap System

The PRS™ Cervical Spine Retention System is a simple and safe solution to effectively position the patient for optimal access and better imaging. Packaged with 4 adhesive pads and 2 (8ft) straps, the Cervical Spine Retention System enhances imaging capabilities of the cervical spine by deflecting and positioning the patient's shoulders down towards the foot of the bed. The large surface area of the pads combined with latex-free, hypoallergenic adhesive reduces skin integrity issues often caused by medical grade tape and alternative adhesives.



Features

- Ideal shoulder deflection for cervical fusion
- Eliminates risk of joint injury caused by pulling wrist restraints for deflection
- System will not release on its own during the procedure
- Adhesive pads are designed to reduce damage to the skin
- Truly one-size-fits-all

Item Number	Description	UOM
CZ-PRS-06	Cervical Spine Retention System 4 Pad, 2 - (96") Strap System	Box of 10

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

22130 R0.1