




TZ Medical
Sparked by your ideas

COMFORT BANDS

Comfort Bands

The TZ Medical Comfort Band offers a soft plastic radial pressure band with ridges on the exterior for fine-tuned adjustments and maximum patient comfort. These bands are available in both sterile and non-sterile varieties, and are also offered in three, color-coded sizes to accommodate all patient types.



Neptune Comfort Bands

The Neptune Comfort Band has all of the great features of the regular Comfort Band, with the addition of a Neptune pad built into the device, to further facilitate hemostasis at the puncture site.



Features

- Bands offer a hands-free approach to hemostasis following needle or sheath removal
- Can be used in cardiac cath labs, radiology departments and dialysis clinics for radial, brachial and hemo-dialysis procedures
- With multiple color-coded lengths to choose from, our bands will fit all body types comfortably
- Cost effective and self-contained hemostasis device



Item #	Description	UOM
Sterile		
CB-8514B	14" blue sterile Comfort Band	Box of 10
CB-8517P	17" purple sterile Comfort Band	Box of 10
Non-Sterile		
CB-7514B	14" blue non-sterile Comfort Band	Each
CB-7517P	17" purple non-sterile Comfort Band	Each
CB-7520G	20" green non-sterile Comfort Band	Each
Neptune		
CB-8514N	14" blue sterile Comfort Band w/ Neptune	Box of 10
CB-8517N	17" purple sterile Comfort Band w/ Neptune	Box of 10

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

10945 R2.0