




TZ Medical
Sparked by your ideas

LEVELUP RADIAL BOARD



The LevelUp Radial Board is designed to elevate the patient's arm so it aligns perfectly with the top of the table mattress, ensuring optimal positioning and comfort during radial catheterization procedures. Its low-profile design easily slides under the table mattress for quick setup and stability. The LevelUp Radial Board can be used on either the left or right side of the table, making it an essential accessory for any cath lab.

- **Extremely durable polypropylene construction**
- **Easy to slide and reposition under the table mattress**
- **Can be used on the right or left side of the table for radial arterial procedures**
- **Improves patient comfort by keeping the arm level with the table mattress**
- **Improves physician posture when accessing the radial sheath, reducing stress and fatigue**



Item #	Description	UOM
LVL-0007	LevelUp Radial Board	1 Each