



TZ Medical
Sparked by your ideas


m·CATH
FLEXY

EXCELLENT CONTROL AND FLEXIBILITY


M-CATH™ Flexy

Microcatheter

- Reinforced tip
- Lengths 135 cm and 150 cm
- Tapered balanced braided shaft
- Durable hydrophilic coating
- Enhanced torsional capacity




Reinforced tip



Great Flexibility



Comfortable Handling




tzmedical.com



800.944.0187



info@tzmedical.com



Technical Information

Product Number	Usable Length	Tip Entry Profile	Dist. Outer Diameter	Prox. Outer Diameter	Guiding Cath. Compatibility	Max. Guidewire Size
M-CATH-03	135 cm	1.65 Fr 0.55 mm	2.46 Fr 0.82 mm	2.79 F 0.93 mm	5 F 1.67 mm	0.014" 0.36 mm
M-CATH-04	150 cm	1.65 Fr 0.55 mm	2.46 Fr 0.82 mm	2.79 F 0.93 mm	5 F 1.67 mm	0.014" 0.36 mm

Manufactured by: Acrostak (Schweiz) AG, Stegackerstrasse 14, CH-8409 Winterthur, Switzerland.
Contact: sales@acrostak.com T: + 41 (0)22 347 68 70 F: + 41 (0)52 233 95 50 www.acrostak.com
M-CATH is trademark of Acrostak (Schweiz) AG

Caution: Federal law restricts this device to sale by or on the order of a physician

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

21732 R0.1