



**TZ Medical**   
ClearBeat Monitoring Service 

# **TRIDENT™** **NANO**

*Holter • Extended Holter*




One Trusted Partner. All of Your Cardiac Monitoring Needs.

# OUR DEVICE

## INNOVATION MEETS SIMPLICITY

The Trident™ Nano is an accurate, reliable, water-resistant, and fully rechargeable monitor designed to capture clear and comprehensive ECG tracings during Holter and Extended Holter studies. 3 channels are recorded from either a lanyard cable or long-wear patch electrode, and its slim, unobtrusive design optimizes patient comfort and compliance.

### FEWER COMPONENTS



### SMARTER MONITORING



Records 3 channels of ECG from a single patch electrode



Fast and simple data retrieval via USB connection



Fully reusable and rechargeable device




3-color LED light indicator for device notifications



Water-resistant for showering

With multiple wear options, 8G of internal storage, and a low-profile design, the Trident™ Nano offers a unique combination of minimal body burden with consistently clear ECG tracings, and the superior diagnostic yield of a 3-channel recording for study durations up to 21 days.

*\*For study durations lasting up to 21 days, multiple patch electrodes are needed.*



MONITORS	ITEM #	DESCRIPTION
Trident™ 30	HPR-0014	Digital Cardiac Monitor for Holter and Extended Holter studies, used with patch electrode and lanyard cable.

One Trusted Partner. All of Your Cardiac Monitoring Needs.

# OUR SERVICE





## Fast & Reliable Arrhythmia Detection

Our certified cardiac technicians provide accurate and reliable analysis of Holter and Extended Holter data 24/7, rapid notification of emergent events, and industry leading turnaround times.

## Customer Support

Our friendly customer service representatives are available to serve you and your patients 24/7. No lengthy teleprompt, no unnecessary hold-times. Just live patient and provider support when you need it, every time.



- *Simplified enrollments*
- *One-step device configuration & reset*
- *Industry leading report turnaround times*
- *Review & e-sign reports from any web-connected device or facility EMR system*

Contact us at [info@clearbeat.com](mailto:info@clearbeat.com) to learn more!

## WHY CHOOSE US

ClearBeat is a comprehensive solution for your practice's Holter, Extended Holter, Cardiac Event, and Mobile Cardiac Telemetry needs. One trusted partner to advance your patients from symptom, to diagnosis, to treatment as quickly, simply, and cost-effectively as possible.

## WHAT SETS US APART

Intuitive, durable, and state-of-the-art monitoring devices coupled with industry-leading report turnaround times, a simple and efficient workflow, and dedication to service after the sale make ClearBeat a clear-cut choice for your ambulatory cardiac monitoring needs.



TZ Medical Inc.  
17750 SW Upper Boones Ferry Rd. #150  
Portland, OR 97224  
Office: 800.944.0187 Fax: 503.639.0239  
www.tzmedical.com  
18874 R1.1

