



**TZ Medical**  
Sparked by your ideas

# PEDIATRIC DEFIB PADS

TZ Medical's Pediatric Defibrillation Electrodes boast the same clinical benefits as our adult pads, however, are designed for your adolescent patients. Our Pediatric Defibrillation Electrodes are constructed to adhere and conform to the patient's body, providing superior current distribution, which results in better clinical outcomes.


Our multi-function electrodes are designed for EP and Cath Labs, Cardiac Surgery, Critical-Care Areas and other related departments. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Compatible with biphasic and monophasic defibrillators, our electrodes offer direct connections to a variety of defibrillators, such as Zoll® and Medtronic® QUIK-COMBO™. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.



## Features

- 60" lead wires for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Unique edge effect technology disperses energy uniformly, reducing burns
- Radiotransparent and radiotranslucent options available

### Radiotransparency



TZ

Competitor

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187




[info@tzmedical.com](mailto:info@tzmedical.com)

19808 R0.6

# Pad Placement Guide



**TZ Medical**  
Sparked by your ideas



## Lead Wire Length

**TZ Sterile P-122-T1: 120"**

**TZ Electrodes: 60"**

**Competitor Electrodes: 42"**

**Competitor Electrodes: 24"**

Medtronic® M1 Connector



Universal T1 Connector



Zoll® Z1 Connector



Product Description	Size	Medtronic® <sup>1</sup>	Universal <sup>2</sup>	Zoll® <sup>3</sup>
Pediatric Radiotranslucent pads	C	P-112-T1 PAD 1111 Adapter*	P-112-T1	P-112-T1 PAD 101PD Adapter*
Pediatric Radiotranslucent pads, sterile	C	P-122-T1 PAD 1111 Adapter*	P-122-T1	P-122-T1 PAD 101PD Adapter*
Pediatric Radiotransparent pads and leadwires	C	P-212-M1		P-212-Z1

**Pad Dimensions:**  
**C=4.25"x2.876"**

**\*Adapters are reusable**

1 – M1 Medtronic models are comparable to Medtronic® Physio-Control QUIK-COMBO™ electrodes.  
 2 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)  
 3 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

**For Pediatric Use:** ≥10kg (22lbs)  
**AED or AED Mode:** ≥25kg (55lbs) ≥8 years old

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

19808 R0.6