



TZ Medical
Sparked by your ideas

EZ HOLDER



The TZ Medical EZ Holder is a versatile device that acts as an anchor point to firmly hold catheters and lead wires in place during procedures. This innovative product will help reduce procedure time as well as improve the efficiency of EP cases and implants.

**Holds your catheters,
so you don't have to!**

Features

- Helps prevent catheter movement during cardioversion, eliminating the need to re-position catheters after therapy is delivered
- Holds the catheter in place for HIS Bundle or any evoked potentials
- Frees up Physician's hands during ablation procedures
- Helps reduce chances of lead dislodgement during implant procedures
- Holds a standard catheter or implant lead that is 5-9 French



Item Number	Description	UOM
-------------	-------------	-----

CT-0010	Catheter and lead wire anchor	Box of 10
---------	-------------------------------	-----------

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

16773 R2.1