



**TZ Medical**  
Sparked by your ideas

# V1-MPC BOOKER BOX



The V1-MPC Patient Companion is a multi-port connector with switching capabilities for connection to defibrillators, EP study equipment, and RF electro-surgical generators. The box provides a stable focal location for multiple wires, cables, and connections.

## The Booker Box System Provides:

- Multi-port connector
- Patient defibrillation
- Cardioversion
- Pacing... and more!



## CALL US TODAY TO SET UP AN EVALUATION



Item Number	Description	UOM
V1-MPC	Booker Box Patient Companion	1 each

## FIND OUT MORE



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

10948 R2.0