




CONTROL THROMBECTOMY FILTER SYSTEM

The Control Thrombectomy Filter is a universal system that is compatible with multiple mechanical thrombectomy devices and provides instant visualization of extracted thrombus without having to leave the sterile field. The transparent housing of the filter is designed with anti-fog material to allow for clear and immediate observation of collected thrombus.



Features

- Provides instant and clear visualization of extracted thrombus on the sterile field
- Filter is compatible with multiple thrombectomy systems
- Fog-resistant chamber traps thrombus for easy visualization

OBSERVE EXTRACTED THROMBUS



WITHOUT LEAVING THE STERILE FIELD

Item Number	Description	UOM
ASP-1039	Control Thrombectomy Filter with on/off switch	Box of 5

FIND OUT MORE

