






TZ Medical
Sparked by your ideas

SHOOT 'N SEE

The TZ Medical Shoot 'N See is a flexible medical grade tape for use during any vascular procedure involving radiography or fluoroscopy. Applied to the skin or drape before the procedure, the radiopaque markings are visible both to the eye and under fluoroscopy, which allows for quick cross-referencing between the inside and outside of a patient's body.

The flexible tape follows patient contours and does not tear under normal use. The Shoot 'N See reduces the need for exposure to contrast media. It also helps with more accurate lesion locating, and more precise stent and catheter sizing.



Features

- Fully visible with or without fluoroscopy
- Arrow indicators to show precise location of lesion(s)
- Allows for quick cross-referencing between inside/outside of patient's body
- Packaged sterile and measures up to 40 cm or 400 mm
- Great for vascular procedures

Item Number	Description	UOM
RCM-0005-10	Shoot 'N See 40 cm tape	Box of 10
RCM-0005-30	Shoot 'N See 40 cm tape	Box of 30
RCM-0005-50	Shoot 'N See 40 cm tape	Box of 50
RCM-0020-10	MM Shoot 'N See 400 mm tape	Box of 10
RCM-0020-30	MM Shoot 'N See 400 mm tape	Box of 30
RCM-0020-50	MM Shoot 'N See 400 mm tape	Box of 50

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

16352 R2.2