



4 Pad, 2 Strap System

The TZ Medical Pannus Retention System is used to safely and effectively position and support excess tissue, allowing easy abdominal access during cesarean sections. Easily applied in one minute, the Pannus Retention System is designed to maintain patient comfort while allowing patient positioning from any angle. The device consists of four adhesive pads and two medical-grade hook & loop straps that can be adjusted after application. Intended for skin applications, the proprietary adhesive in the pads is residue-free and minimizes the risk of skin tears and degradation typically seen with the use of tape and competing products.



4 - 5.5" x 7.5" Pads



2 - 48" Strap

Features

- Hypoallergenic, latex & residue free
- Adhesive pads will not tear or damage skin
- Ideal for High-BMI Cesarean Sections and Trans Aortic Valve Replacements (TAVR)
- Completely adjustable system - will fit any sized patient
- Can be applied in less than one minute
- Easy to reposition after application
- Sterile and designed for surgical procedures



Safely retains large amounts of tissue



Item Number	Description	UOM
CZ-PRS-04	Pannus Retention System 4 pad, 2 strap system, sterile	10/Box

FIND OUT MORE

