



TZ Medical
Sparked by your ideas

PANNI™ POSTPARTUM RETRACTOR

The Panni™ Postpartum Retractor is a postoperative pannus management system designed to retract and position excess tissue away from the incision site after C-section procedures. The Panni™ can be easily applied by one person and is completely adjustable to ensure optimal patient comfort. This device helps to facilitate wound health by exposing the incision site to air and can be worn for up to 3 days while laying, standing, and showering.



Features

- Easily applied by one person
- Shower compatible
- Can be worn for up to 3 days
- Completely adjustable
- Can be applied in seconds
- Residue-free
- Saves time



Water Resistant



Readjustable

| Item Number | Description | UOM |
|-------------|-------------|-----|
|-------------|-------------|-----|

| | | |
|--------|----------------------|-----------|
| CZ-PPR | Postpartum Retractor | Box of 10 |
|--------|----------------------|-----------|

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

19566 R0.2