




Thread Comparison




TZ Medical
Titanium
Distraction Pin


Stainless Steel
Distraction Pin

Better ultimate strength means better performance and safety.

Ultimate Strength Comparison



Item Number	Description	UOM
DP-12-TB	12mm Titanium Distraction Pin Blue	Box of 10
DP-14-TY	14mm Titanium Distraction Pin Yellow	Box of 10
DP-16-TG	16mm Titanium Distraction Pin Green	Box of 10
DP-18-TP	18mm Titanium Distraction Pin Purple	Box of 10



Distraction Pin Accessories

DP-S	Caspar Screwdriver Universal
DP-72-CDR	Caspar Distractor Right
DP-74-CDL	Caspar Distractor Left

Features

- Triple Fluted design for easier insertion into bone
- Wider diameter than typical distraction pins
- Color-coded sizes for faster identification
- Made with 100% titanium
- Cost effective
- Safer and sturdier than pins made from stainless steel
- Compatible with all current ACF systems

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

11946 R2.1