

COMFORT HALO



TZ Medical
Sparked by your ideas


The TZ Medical Comfort Halo is designed to enhance safety and comfort for patients undergoing conscious sedation.

The polycarbonate base is compatible with all bed types, allowing it to easily slide underneath the imaging table mattress pad. The steel and vinyl tubing can then be manipulated into a multitude of different draping positions, thereby keeping the drape out of the patient's face and allowing complete visualization of the patient at all times.



For All Procedures Using Conscious Sedation Including:

- PPM/ICD/Bi-V ICD implants
- Chemo/Port placements
- CVP and PICC line insertion
- Antegrade femoral angiography
- Rt Heart Caths from Subclavian Approach



Features

- Increases patient comfort by keeping anesthesia drape off of patient's face during conscious sedation
- Ideal for all departments using radiographic imaging and sterile technique
- Offers an unobstructed view of patient during procedure
- The flexible design may be molded to the desired shape or height
- Visual reminder to perform "Time Out", enhancing patient safety/ improving staff compliance with JCAHO requirements
- Compatible with any bed type, no arm rails necessary

Item Number	Description	UOM
CZ-009	Comfort Halo	1 each

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

12596 R1.0