




The Patient Immobilizer allows the patient to be secured during imaging procedures in order to limit excess movement, thereby eliminating motion artifact. This greatly reduces the need for repeat subtraction runs, which in turn decreases patient and staff radiation exposure.

The patient's wrists can be safely attached to the leg wrap by using the quick snap and rapid adjustment webbing, which can also be quickly undone as soon as the procedure is complete. For additional safety, a second wrap may also be secured around the patient's chest.



Features

- Improves safety for patient and staff during imaging procedures
- Eliminates motion artifact caused by patient movement
- Cuts down on the need for repeat subtraction runs.
- Quick snap and rapid adjustment webbing for fast immobilization of patient
- Patient may be effectively secured during procedures where conscious sedation is used

The Patient Immobilizer XXL Extender allows patients of any size to be safely secured during imaging procedures. The XXL Extender attaches easily to the Patient Immobilizer and offers a much greater range of size adjustment.



Item Number	Description	UOM
CZ-106	Patient Immobilizer	1 kit
CZ-107	GE Patient Immobilizer	1 kit
CZ-008	Disposable Wrist Wrap	Box of 10
CZ-012	XXL Extender	1 each

FIND OUT MORE

