




POLYPROPYLENE COMFORT CRADLE

Comfort Cradle Dimensions



The Polypropylene Comfort Cradle is a versatile product that aids both the patient and the operator during procedures. The cradle's extended width improves comfort, making long procedures more tolerable for the patient.

The Polypropylene Comfort Cradle also offers the ability to position the Flat Detector (or II) in direct contact with the patient in all views, reducing secondary and scatter radiation. The result is superior image quality and reduced radiation exposure for the staff.

The dimensions of the cradle provide superior support and comfort during vascular procedures when compared with smaller and less robust arm boards.

Features

- Patient is more accessible to the physician
- Extra wide dimensions for added comfort
- Reduces secondary and scatter radiation
- Extremely durable Polypropylene construction
- Superior image quality
- Comes with a six-month warranty



Item Number	Description	UOM
CZ-100-PP	Set of 2 Polypropylene Comfort Cradle arm boards	1 kit

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

17905 R1.0