



TZ MEDICAL StarTrol LED MEDICAL LIGHTS

TZ MEDICAL STARTROL MEDICAL LIGHTS

StarTrol LED Medical Lights produce a bright white light, (94 CRI and 4700K). They are ideal for true color rendition, and are cool and comfortable to work under, with virtually no radiant heat being generated by the LED's. The lights are energy efficient, using less than one-third the energy of a halogen light, and require minimal maintenance. StarTrol LED Medical Lights are durable, precision machined in the USA and feature a 10-year warranty.

TZ Medical offers a wide variety of light configurations, including ceiling mounted, wall mounted, and mobile rolling units, which are sure to fit into any work space.

CEILING MOUNTED UNITS



*StarTrol 8X4-GS1
Ceiling Mounted Light
Item# 8X4-GS1*



*StarTrol 4X4-GS1
Ceiling Mounted Light
Item# 4X4-GS1*

DUAL ARM CEILING MOUNTED UNITS



*StarTrol 8X4-GS2
Ceiling Mounted Light
Item# 8X4-GS2*



*StarTrol 4X4-GS2
Ceiling Mounted Light
Item# 4X4-GS2*

PRODUCT FEATURES

- Energy Efficient & Eco-friendly
- Bright, natural white light
- 10-year limited warranty
- Cool to work under
- 40,000 hours of LED life
- Infinite variable positioning
- Multi-level dimming



TZ MEDICAL StarTrol LED MEDICAL LIGHTS

MOBILE UNITS



*StarTrol 4X4-GSM
Mobile Examination Light
Item# 4X4-GSM*



*StarTrol 8X4-GSM
Mobile Examination Light
Item# 8x4-GSM*

WALL MOUNTED UNITS



*StarTrol 8X4-GSW
Wall Examination Light
Item# 8X4-GSW*



*StarTrol 4X4-GSW
Wall Examination Light
Item# 4X4-GSW*

*More lighting configurations are available.
For a complete list of our medical lighting,
and to check out our other product lines,
visit us online at www.tzmedical.com or call
us at 800.944.0187*

17750 SW Upper Boones Ferry Rd. #150
Portland, OR 97224
Office: 800.944.0187 Fax: 503.639.0239
www.tzmedical.com



TZ Medical
Sparked by your ideas