



TZ Medical
Sparked by your ideas

TRUBRITE



Adhesive bottom for easy placement on various surfaces



Simple pull tab activates the light



Features

- Sterile, compact & disposable design that weighs less than 4 grams
- Provides continuous light for over 2 hours
- Adjustable LED light, providing optimal visibility for any procedure
- Self-adhesive for effortless mounting on multiple surfaces
- Ideal for pacemaker pockets with Army-Navy retractors



LED light can be positioned straight or bent for optimal illumination

Item Number	Description	UOM
ARL-1010	Adhesive Retractor Light	Box of 10

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

18618 R1.2