



# FLAT RADIAL BOARD



*The dimensions of the Flat Radial Board are 44.4" x 22.9", offering plenty of space and positioning options.*

The Flat Radial Board is a lightweight carbon fiber board specifically shaped to allow convenient radial access in conjunction with medical imaging and patient comfort. By combining radial access with a landing zone for equipment, the Flat Radial Board takes radial access to a whole new level.

The Flat Radial Board can be used for radial access on either the right or left side.

*The Flat Radial Board slides under the mattress for increased stability.*



## Features

- Can be used on either the right or left side for radial access
- Large landing area for equipment
- Lightweight carbon fiber construction (No foam core)
- Invisible under fluoroscopy
- Can be used for all diagnostic and interventional imaging procedures

Item Number	Description	UOM
CZ-400	Flat Radial Board	1 each

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

15508 R1.0