



The Flipper is a versatile radiation shield designed to keep physicians safe without getting in the way of the procedure. The Flipper's 19"x16" versatile design allows for use on either the right or left side for femoral groin access or subclavian access. Convenient adhesive strips on either side of the shield allow it to be secured in place to avoid slipping.

Lead free, and attenuating up to 75% of radiation at 90 kVp, the Flipper is ideal for Interventional Radiology, Cardiac Cath Labs, Electrophysiology and more.



## Features

- Attenuates up to 75% of radiation at 90 kVp
- Reduces scatter radiation exposure to operators and staff
- Can be used for femoral groin access or subclavian access
- Adhesive on both sides to prevent slipping
- Lead equivalence 0.125mm
- Perfect for Interventional Radiology, Cardiac Cath, Electrophysiology & more

Item Number	Description	UOM
CZ-007	Flipper Radiation Shield	Box of 15

**FIND OUT MORE**



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

15104 R2.2