




The TZ Medical URPro is a universal system that directly replaces your current MAVIG™ shields with high-quality and long-lasting radiation protection panels. The URPro panels are uniquely designed to match your current MAVIG™ frame system\*, allowing you to replace any and/or all damaged panels through a simple installation process. More importantly, each panel is well-structured with fabric binding around the edges to reduce common issues with fraying and deterioration caused by long-term usage.

\*TZ Medical only supplies URPro Radiation Protection Shields/Panels. All other components are manufactured and trademarked products of the MAVIG™ corporation.

## Features

- Attenuates 98% of secondary and scatter radiation at 80 kVp
- Lead equivalent of 0.50 mm
- Enhanced lower body radiation protection
- Directly Compatible with current MAVIG™ bedside drape framing systems
- Simple installation process
- Can replace any number of damaged panels
- Fabric binding and structured frame reduces long-term fraying and deterioration
- Long-lasting solution at a cost-effective price



**OLD PANELS**

**NEW PANELS**

## FIND OUT MORE



[tzmedical.com](http://tzmedical.com)



800.944.0187



[info@tzmedical.com](mailto:info@tzmedical.com)

20583 R0.2

Item Number	Description	UOM
BDR-1002	URPro Replacement Panel	1 panel
BDR-1002-5	URPro Replacement Set	5 panels