



TZ Medical
Sparked by your ideas


STERILE DEFIB PADS



Our multi-function sterile electrodes are ideal for a variety of procedures, including TAVR, AVR, Mini-Maze, ICD Implants, Laser Lead Extraction and much more. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Compatible with biphasic defibrillators, our sterile electrodes offer direct connection to Zoll® defibrillators. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.



Radiotransparency



TZ

Competitor

Features

- Lead wire lengths available in 60" and 120" for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Radiotransparent and radiotranslucent options available

FIND OUT MORE



tzmedical.com



800.944.0187




info@tzmedical.com

13622 R2.9

Pad Placement Guide



TZ Medical
Sparked by your ideas



Lead Wire Length

TZ Sterile P-121-T1, P-122-T1: 120"

TZ Sterile P-224-Z1, P-224-T1: 60"

Competitor Electrodes: 24"

Universal T1
Connector



Zoll® Z1
Connector



Product Description	Size	Medtronic®	Universal ¹	Zoll® ²
Adult Radiotranslucent, sterile	A	P-121-T1 PAD 1111 Adapter*	P-121-T1	P-121-T1 PAD 101PD Adapter*
Adult Radiotransparent pads and leadwires, large posterior pad, sterile	A/B	P-224-T1 PAD 1111 Adapter*	P-224-T1	P-224-Z1
Pediatric Radiotranslucent pads, sterile	C	P-122-T1 PAD 1111 Adapter*	P-122-T1	P-122-T1 PAD 101PD Adapter*

Pad Dimensions:

A=6"x 4.25" B=6"x9" C=4.25"x2.876"

*Adapters are reusable

1 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)

2 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

For Pediatric Use: ≥25kg (55lbs)

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

13622 R2.9