



TZ Medical
Sparked by your ideas

Medical Mattress Catalog





Contents

-  Stretcher & Gurney 5-20
-  Operating Tables 21-24
-  X-Ray, Angio & Cath Lab . . . 25-26
-  Overlay/Comfort Pads 27-30
-  General Patient/ICU 31-35
-  Bariatric 36-38
-  V36 Deck Bed 39
-  Side Rail Pads 40-41
-  Labor and Delivery 42-46
-  Pediatric 47-50

All of our Products are Made in the U.S.A.

Our extensive range includes hospital mattresses, stretcher pads, overlays, and disaster relief mattresses.

Our American-made mattresses are created for comfort, support, durability, quality and value.



INVEST IN A QUALITY MATTRESS

Industry Leading Comfort

Our pressure mapped custom foam layering ensures a level of comfort to provide a restful night's sleep.

Durability You Can Trust

Advanced stability layers provide long lasting comfort without sagging, ensuring a consistent, hard-working performance every time it's used.

Easily Cleaned Removable Covers

Our covers that are expertly made to provide a safe, comfortable, and durable solution for the healthcare industry.

Quality Assurance

Proudly American-made and veteran owned, we stand behind our products with a 12-month warranty.

Replacing Your Mattresses

By staying proactive in maintaining your hospital mattresses, you are safeguarding the well-being of your patients and promoting a hygienic environment within your healthcare facility. Follow these tips to ensure the integrity of your mattress.



Develop an Inspection Plan

Create an inspection plan for all hospital bed mattresses and mattress covers in your facility.

Follow the manufacturer's guidelines for expected lifespan and any additional recommendations. Replace when necessary.

Contact the mattress cover manufacturer for any additional questions not covered here.



Inspect

Regularly check each hospital mattress cover for any visible signs of damage or wear such as cuts, tears, cracks, pinholes, snags, or stains.

Routinely remove the hospital bed mattress cover and check inside surface. Once the mattress cover is removed, inspect the mattress for wet spots, staining, or signs of damage or wear. Check all sides and the bottom of the mattress.

Be aware that it may be difficult to identify damaged or soiled mattresses without removing the mattress cover first. Mattress covers tend to be dark in color, making it hard to see what lies underneath.



Remove and Replace

Remove any damaged, worn, or visibly stained hospital bed mattress according to the healthcare facility's procedures and manufacturer's instructions.

Immediately replace any mattress cover with visible signs of stains, damage or wear to reduce the risk of infection to patients.



Maintain

Clean and disinfect undamaged mattress covers according to the manufacturer's guidelines.

DO NOT stick needles into a hospital bed mattress through the mattress cover.



Cleaning and Disinfecting

GENERAL CLEANING:

- Soils or Stains: Use neutral soap suds and lukewarm water. Do not use harsh cleansers, solvents or detergents.
- Hard-to-Clean Spots: Use standard household/vinyl cleansers, a soft bristle brush and gentle agitation on troublesome spots or stains. Pre-soak heavy, dried-on soil.

DISINFECTING:

- Phenolic disinfectants, such as “Simple Green d PRO 5”, are the best choice for vinyl products, though properly diluted quaternaries are also acceptable. Phenolic should be avoided on urethane.
- Ammonia disinfectants are recommended for urethane products. Quaternary/Isopropyl disinfectants are recommended for urethanes, but not recommended for vinyl.
- All NAMC fabrics are cleaned with a 1:30 dilution of “Simple Green d PRO 5” There is no harmful effect on the fabric when used in this dilution ratio.

SPECIFIC COMMERCIAL DISINFECTANTS:

- There are many popular disinfectants being used in care facilities. Cavi-Wipes and Virox5 are two of these that have been tested by Herculite and found to be acceptable for use with NAMC healthcare fabrics when used according to manufacturer’s recommendations.
- Other specific products can be tested by NAMC when the request is accompanied by an unopened product sample. Contact your Sales Representative to coordinate this. Results are communicated in approximately 4-6 weeks.

LAUNDERING:

- Machine washing or drying vinyl-laminated or rubber-coated NAMC fabrics is not recommended. Laundering may substantially decrease the useful life of the fabric.
- Polyurethane-coated Sure-Chek fabrics: machine-wash with detergent up to 120°. (Some surface wrinkling may occur. This wrinkling has no adverse effect on the fabric properties). Hang or tumble dry thoroughly at 140° before storage. Using Hypochlorite as recommended by the Center for Disease Control, there is no harmful effect on the fabric when used in this dilution ratio.

PRECAUTIONS:

- Disinfectants applied at full concentration or in highly concentrated solutions may damage the fabric and will decrease the useful life of Sure-Chek fabrics.
- Iodophor type disinfectants used on NAMC fabrics may result in staining. If stains from Iodopher (i.e. iodine bearing) type disinfectants are not treated with a diluted (10:1) bleach solution within 20 minutes of application or spillage, staining will occur.
- Storage of fabrics that have not been dried thoroughly before folding and storing can promote the growth of mildew.

LIST OF CHEMICALS TESTED BY NAMC:

Believed to be safe to use on NAMC brand healthcare fabrics when used according to manufacturer’s directions.

- Simple Green d PRO 5 (used by NAMC)
- Virox 5 Cleaner
- Oxivir® TB Disinfectant Cleaner
- CaviWipes Disinfecting Towels
- ProCure V - Active Ingredient: Chlorine Dioxide
- PDI Sani-Cloth HB Germicidal Disposable Wipe
- Clorox Healthcare Bleach Germicidal Wipes
- Super Sani-Cloth
- Pro-Q256 Cleaners

Stretchers



Fabrics

Bed and stretcher mattresses come in a variety of materials. Our Ultra Comfort stretchers feature Nymed 6™ healthcare grade nylon top covers and a non-slip, gripable, four way stretch, bottom covers. Standard stretcher models feature Vynmed™ healthcare vinyl top and bottom covers.

	Nymed™ Nylon	Vynmed™ Vinyl
Breathable	●	●
Keeps Patients Cool and Dry	●	●
Anti-bacterial, Anti-fungal, Hypo allergenic	●	●
Non-Conductive, Antistatic	●	●
Stain Resistant	●	●
Fluid Resistant	●	●
Non-allergenic, Phthalate Free, Latex Free, Non-irritating	●	●
Odor Resistant	●	●
Easy to Clean	●	●
Our overall most durable fabric	●	●



Color Identifier

Color Identifier (tagging) at the mattress foot is now included as a standard feature on all Stretcher Mattresses.

Distinguish old mattresses from new or use for departmental identification. Standard yellow identifier (tagging) included, optional colors (red or blue) available upon request.

Choosing the Right Mattress

Stretcher & Gurney Mattresses feature multi-layered, CertiPUR-US foam. The removable covers are fluid proof, anti-bacterial, anti-fungal, stain resistant and come with a protective fluid flap over the heavy-duty zipper located under the mattress (the zipper is located at the foot) to ensure no fluid contamination. All mattresses include a standard yellow color identifier at foot.

Features	Standard	Ultra-Comfort	Bariatric
Cool Gel Infused Memory Foam		●	
Multi-layer High-resiliency Foam	●	●	●
Water Resistant Cover	●	●	●
Secondary Fluid Barrier (protects foam if outer layer is compromised)	●	●	●
Removable Cover	●	●	●
Fluid Protection Flap Over Zipper	●	●	●
Color Identifier (tagging)	●	●	●
Easy to Clean	●	●	●
Shear Liner (optional feature)	●	●	●
Latex Free	●	●	●
FirmEdge™ Passive Restraint System		●	●
Weight Capacity	300lbs	300lbs	500lbs
Mattress Weight	20lbs	25lbs	40lbs
Warranty	12/mo.	12/mo.	12/mo.
Made in the USA	●	●	●
CertiPUR-US Certified Foam	●	●	●

Standard Mattress Features



Heavy-duty Zipper

On bottom allows for internal mattress inspection & easy replacement.



Fluid Protection Flap Over Zipper

Designed to direct fluids over and away from the zipper area.



Flame Retardant

Meets California Technical Bulletin 121 and CFR 1632.



Secondary Fluid Barrier

Assembled between foam and cover providing unique two-way protection.

- Protects foam if outer cover is compromised
- Protects cover by preventing foam's plasticizers from embrittling the fabric, extending cover life



CertiPUR-US Certified Foam

- Made without ozone depleters, PBDEs, TDCPP or TCEP ("Tris") flame retardants
- Made without mercury, lead and other heavy metals
- Made without formaldehyde
- Made without phthalates regulated by the Consumer Product Safety Commission
- Low VOC (Volatile Organic Compound) emissions for indoor air quality (less than 0.5 parts per million)



Color Identifier (tagging)

Standard yellow identifier (tagging) included, optional colors (red or blue) upon request.



12 Month Standard Warranty

Standard warranty covers material and workmanship. Extended warranty available.



Made in the USA - Veteran Owned

Quality workmanship using locally sourced materials.

Unique Features



Removable Nymed™ Cover

Flame resistant, water resistant, anti-bacterial, anti-fungal, antistatic, non-allergenic, non-irritating, and is latex free.



Removable Vynmed™ Cover

Durable, breathable Vynmed™ healthcare vinyl. Soft-supple feel offers exceptional comfort. Anti-bacterial, anti-fungal, water resistant, stain resistant, antistatic, non-allergenic, non-irritating and is latex free.



FirmEdge™ Passive Restraint System

Provides support on both sides of the mattress during entry and egress. (On ultra comfort model only)



Patient Weight Limit

Depending on the model selected our mattress provides support for patients up to 300lbs. for the standard and ultra comfort.



Bariatric model supports up to 500lbs. and 1000lbs.



Cool Gel Memory Foam

Plush and luxurious, evenly distributes body weight, infused liquid gel for a cooler sleep experience, wicks heat away from body. (On ultra comfort model only)



Multi-Layer Performance Foam

Soft surface/firm base for comfort and support.




Standard Stretcher

Multi-layered, Medical Grade,
High Resiliency, Urethane Foam

Our economical, standard stretcher replacement mattresses meet or exceed OEM specifications.

4" thick CertiPUR-US foam provides comfortable support for individuals up to 300lbs.

- Soft conforming layer
- Supportive firm foam base layer



Color Identifier located at foot end.

Distinguish old mattresses from new or use for departmental identification. Standard yellow identifier (tagging) included, optional colors available upon request.

Vynmed™ Healthcare Vinyl Cover

- Durable, soft-supple feel offers exceptional comfort
- Anti-bacterial, anti-fungal, water resistant, stain resistant
- Non-allergenic, non-irritating, conductive
- Phthalate free, latex free and odor resistant


Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area

NEW! Non Slip Bottom Material



- Multi-layer Performance Foam
- 300 Patient Weight Limit
- Removable Vynmed™ Cover
- Heavy-duty Zipper
- Fluid Protection Flap Over Zipper
- Secondary Fluid Barrier
- Universal Velcro Location For Easy Installation
- Color Identifier
- CertiPUR-US Foam
- 12 Month Warranty
- Made in the USA



Ultra-Comfort Stretchers

Stretcher & Gurney Mattresses

Cool Gel Memory Foam, Non Slip Bottom

Nymed 6™ stretchable nylon cover accompanied with our multi-layered, ultra comfort, medical grade, high resiliency, urethane foam.

4" or 5" thick CertiPUR-US foam provides comfortable support for individuals up to 300lbs.

- Cool Gel infused memory foam top layer
- Performance soft foam internal layer
- Supportive firm foam base layer

Color Identifier located at foot end.

Distinguish old mattresses from new or use for departmental identification. Standard yellow identifier (tagging) included, optional colors available upon request.

FirmEdge™ Passive Restraint System

- Provides support on both sides of the mattress during entry and egress

Nymed 6™ Healthcare Grade Cast Urethane Top Cover

- Water resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Non Slip Bottom Material


- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Phthalate free, latex free and odor resistant
- Keeps mattress on stretcher

Secondary Fluid Barrier



- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area




NEW! Non Slip Bottom Material



Cool Gel Memory Foam




FirmEdge™ Passive Restraint




300 Patient Weight Limit




Removable Cast Urethane Cover




Heavy-duty Zipper




Fluid Protection Flap Over Zipper




Secondary Fluid Barrier




Non Slip Bottom




Color Identifier



CertiPUR-US Foam










12 Month Warranty










Made in the USA








Stryker Stretchers

Stryker	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (26" Width)	Transport	Stryker 700	TZ-700-4	TZ-700-4-UC	TZ-700-5-UC
 (24" and 26" Widths)	Transport	Stryker 721	TZ-721-24-4 TZ-721-26-4	TZ-721-24-4-UC TZ-721-26-4-UC	TZ-721-24-5-UC TZ-721-26-5-UC
 (26" and 30" Widths)	Transport	Stryker 737	TZ-737-26-4 TZ-737-30-4	TZ-737-26-4-UC TZ-737-30-4-UC	TZ-737-26-5-UC TZ-737-30-5-UC
 (26" and 30" Widths)	Transport	Stryker 738	TZ-738-26-4 TZ-738-30-4	TZ-738-26-4-UC TZ-738-30-4-UC	TZ-738-26-5-UC TZ-738-30-5-UC
 (24" Width)	Gurney	Stryker 916	TZ-916-4	TZ-916-4-UC	TZ-916-5-UC
 (26" Width)	Transport	Stryker 919	TZ-919-4	TZ-919-4-UC	TZ-919-5-UC
 (24" and 26" Widths)	ED II	Stryker 926	TZ-926-24-4 TZ-926-26-4	TZ-926-24-4-UC TZ-926-26-4-UC	TZ-926-24-5-UC TZ-926-26-5-UC
 (26" and 29" Widths)	PACU	Stryker 946	TZ-946-26-4 TZ-946-29-4	TZ-946-26-4-UC TZ-946-29-4-UC	TZ-946-26-5-UC TZ-946-29-5-UC
 (24" Width)	Emergency	Stryker 952	TZ-952-4	TZ-952-4-UC	TZ-952-5-UC








Stryker Stretchers CONT.

Stryker	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (24" Width)	Emergency	Stryker 960	TZ-960-4	TZ-960-4-UC	TZ-960-5-UC
 (24" Width)	Emergency	Stryker 966	TZ-966-4	TZ-966-4-UC	TZ-966-5-UC
 (26" and 30" Widths)	VIP Table	Stryker 974	TZ-974-26-4 TZ-973-30-4	TZ-974-26-4-UC TZ-974-30-4-UC	TZ-974-26-5-UC TZ-974-30-5-UC
 (24" Width)	Advantage	Stryker 1001	TZ-1001-4	TZ-1001-4-UC	TZ-1001-5-UC
 (26" and 30" Widths)	Trauma	Stryker 1002	TZ-1002-26-4 TZ-1002-30-4	TZ-1002-26-4-UC TZ-1002-30-4-UC	TZ-1002-26-5-UC TZ-1002-30-5-UC
 (26" and 30" Widths)	ER PACU Transport	Stryker 1007	TZ-1007-26-4 TZ-1007-30-4	TZ-1007-26-4-UC TZ-1007-30-4-UC	TZ-1007-26-5-UC TZ-1007-30-5-UC
 (26" Width)	Advantage	Stryker 1009	TZ-1009-4	TZ-1009-4-UC	TZ-1009-5-UC
 (26" Width)	ER & Transport	Stryker 1010	TZ-1010-4	TZ-1010-4-UC	TZ-1010-5-UC
 (26" and 30" Widths)	Transport Big Wheel	Stryker 1015	TZ-1015-26-4 TZ-1015-30-4	TZ-1015-26-4-UC TZ-1015-30-4-UC	TZ-1015-26-5-UC TZ-1015-30-5-UC

Stryker Stretchers CONT.

Stryker	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (26" Width)	Trauma	Stryker 1020	TZ-1020-4	TZ-1020-4-UC	TZ-1020-5-UC
 (26" and 30" Widths)	Prime Zoom M Series	Stryker 1025	TZ-1025-26-4 TZ-1025-30-4	TZ-1025-26-4-UC TZ-1025-30-4-UC	TZ-1025-26-5-UC TZ-1025-30-5-UC
 (26" Width)	L & D	Stryker 1060 (2 Piece Set)	TZ-1060-4	TZ-1060-4-UC	TZ-1060-5-UC
 (26" Width)	L & D	Stryker 1061 (2 Piece Set)	TZ-1061-4	TZ-1061-4-UC	TZ-1061-5-UC
 (24" Width)	EYE	Stryker 1067	TZ-1067-4	TZ-1067-4-UC	TZ-1067-5-UC
 (26" Width)	EYE	Stryker 1068 (Body Section)	TZ-1068-4	TZ-1068-4-UC	TZ-1068-5-UC
	EYE	Stryker 1068 (Head Section - Flat)	TZ-1068-4-H-F	TZ-1068-4-H-F-UC	TZ-1068-5-H-F-UC
	EYE	Stryker 1068 (Head Section - Concave)	TZ-1068-4-H-C	TZ-1068-4-H-C-UC	TZ-1068-5-H-C-UC





Stryker Stretchers CONT.

Stryker	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (26" Width)	EYE	Stryker 1069 (Body Section)	TZ-1069-4	TZ-1069-4-UC	TZ-1069-5-UC
	EYE	Stryker 1069 (Head Section - Flat)	TZ-1069-4-H-F	TZ-1069-4-H-F-UC	TZ-1069-5-H-F-UC
	EYE	Stryker 1069 (Head Section - Concave)	TZ-1069-4-H-C	TZ-1069-4-H-C-UC	TZ-1069-5-H-C-UC
 (26" Width)	Fluoroscopy	Stryker 1080	TZ-1080-4	TZ-1080-4-UC	TZ-1080-5-UC
 (26" Width)	EYE	Stryker 1089 (Body Section)	TZ-1089-4	TZ-1089-4-UC	TZ-1089-5-UC
	EYE	Stryker 1089 (Head Section - Flat)	TZ-1089-4-H-F	TZ-1089-4-H-F-UC	TZ-1089-5-H-F-UC
	EYE	Stryker 1089 (Head Section - Concave)	TZ-1089-4-H-C	TZ-1089-4-H-C-UC	TZ-1089-5-H-C-UC




Stryker Stretchers CONT.

Stryker	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (26" and 30" Widths)	5th Wheel Prime	Stryker 1105	TZ-1105-26-4 TZ-1105-30-4	TZ-1105-26-4-UC TZ-1105-30-4-UC	TZ-1105-26-5-UC TZ-1105-30-5-UC
 (26" and 30" Widths)	Prime Series Big Wheel	Stryker 1115	TZ-1115-26-4 TZ-1115-30-4	TZ-1115-26-4-UC TZ-1115-30-4-UC	TZ-1115-26-5-UC TZ-1115-30-5-UC
 (26" and 30" Widths)	Zoom Prime	Stryker 1125	TZ-1125-26-4 TZ-1125-30-4	TZ-1125-26-4-UC TZ-1125-30-4-UC	TZ-1125-26-5-UC TZ-1125-30-5-UC
 (26" and 29" Widths)	Renaissance	Stryker 1210	TZ-1210-26-4 TZ-1210-29-4	TZ-1210-26-4-UC TZ-1210-29-4-UC	TZ-1210-26-5-UC TZ-1210-29-5-UC
 (26" and 29" Widths)	Renaissance	Stryker 1211	TZ-1211-26-4 TZ-1211-29-4	TZ-1211-26-4-UC TZ-1211-29-4-UC	TZ-1211-26-5-UC TZ-1211-29-5-UC
 (29" Width)	Advantage	Stryker 1501	TZ-1501-4	TZ-1501-4-UC	TZ-1501-5-UC
 (26" and 29" Widths)	PACU	Stryker 1509	TZ-1509-26-4 TZ-1509-29-4	TZ-1509-26-4-UC TZ-1509-29-4-UC	TZ-1509-26-5-UC TZ-1509-29-5-UC
 (26" Width)	Advantage	Stryker 1510	TZ-1510-4	TZ-1510-4-UC	TZ-1510-5-UC






Stryker Stretchers CONT.

Stryker	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (26" and 30" Widths)	Advantage	Stryker 1550	TZ-1550-26-4 TZ-1550-30-4	TZ-1550-26-4-UC TZ-1550-30-4-UC	TZ-1550-26-5-UC TZ-1550-30-5-UC
 (26" Width)	Advantage	Stryker 1710	TZ-1710-4	TZ-1710-4-UC	TZ-1710-5-UC
 (26" and 29" Widths)	Advantage	Stryker 1711	TZ-1711-26-4 TZ-1711-29-4	TZ-1711-26-4-UC TZ-1711-29-4-UC	TZ-1711-26-5-UC TZ-1711-29-5-UC
 (24" Width)	Eye Stretcher	Stryker 5051	TZ-5051-4	TZ-5051-4-UC	TZ-5051-5-UC

Pedigo Stretchers

Pedigo	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
	Transport	Pedigo 511/540/550	TZ-5829002	TZ-5828002-UC	TZ-5829002-UC
	Transport	Pedigo 516/545/555	TZ-5832002	TZ-5831002-UC	TZ-5832002-UC
	Transport	Pedigo 530	TZ-5803002	TZ-5802002-UC	TZ-5803002-UC




Midmark Stretchers

Midmark	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
	Transport	Midmark 511	TZ-511-4	TZ-511-4-UC	TZ-511-5-UC
	Transport	Midmark 516	TZ-516-4	TZ-516-4-UC	TZ-516-5-UC
	Universal	Midmark 540	TZ-540-4	TZ-540-4-UC	TZ-540-5-UC
	PACU	Midmark 545	TZ-545-4	TZ-545-4-UC	TZ-545-5-UC
	Universal	Midmark 550	TZ-550-4	TZ-550-4-UC	TZ-550-5-UC





Hill-Rom Stretchers

Hill-Rom	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
	GPS	Hill-Rom 880	TZ-880-4	TZ-880-4-UC	TZ-880-5-UC
	GPS	Hill-Rom 881	TZ-881-4	TZ-881-4-UC	TZ-881-5-UC
	GPS	Hill-Rom 882	TZ-882-4	TZ-882-4-UC	TZ-882-5-UC
	GPS	Hill-Rom 883	TZ-883-4	TZ-883-4-UC	TZ-883-5-UC
	GPS	Hill-Rom 885	TZ-885-31-4	TZ-885-31-4-UC	TZ-885-31-5-UC
	GPS	Hill-Rom 886	TZ-886-31-4	TZ-886-31-4-UC	TZ-886-31-5-UC
 (26" and 30" Widths)	TranStar	Hill-Rom 8000	TZ-8000-26-4 TZ-8000-30-4	TZ-8000-26-4-UC TZ-8000-30-4-UC	TZ-8000-26-5-UC TZ-8000-30-5-UC
	DuraStar	Hill-Rom 8005	TZ-8005-4	TZ-8005-4-UC	TZ-8005-5-UC
	TranStar Surgical	Hill-Rom 8010	TZ-8010-4	TZ-8010-4-UC	TZ-8010-5-UC
	TranStar Electric	Hill-Rom 8020	TZ-8020-4	TZ-8020-4-UC	TZ-8020-5-UC

Hill-Rom Stretchers CONT.

Hill-Rom	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
	TranStar Fixed	Hill-Rom 8030	TZ-8030-4	TZ-8030-4-UC	TZ-8030-5-UC
	TranStar Trauma	Hill-Rom 8040	TZ-8040-4	TZ-8040-4-UC	TZ-8040-5-UC
	TranStar OB/Gyn	Hill-Rom 8050	TZ-8050-4	TZ-8050-4-UC	TZ-8050-5-UC







Hausted Stretchers

Hausted	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
 (28" and 30" Widths)	Horizon Series Ameriglide & Retracto	Hausted 428	TZ-428-28-4 TZ-428-30-4	TZ-428-28-4-UC TZ-428-30-4-UC	TZ-428-28-5-UC TZ-428-30-5-UC
	Horizon Youth Series	Hausted 462	TZ-462-4	TZ-462-4-UC	TZ-462-5-UC
	Converge II	Hausted 471	TZ-471-4	TZ-471-4-UC	TZ-471-5-UC
	Converge II	Hausted 472	TZ-472-4	TZ-472-4-UC	TZ-472-5-UC

Hausted Stretchers CONT.

Hausted	Stretcher	Model	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
	Surgi-Stretcher	Hausted 575	TZ-575-4	TZ-575-4-UC	TZ-575-5-UC
	Transport	Hausted 600	TZ-600-4	TZ-600-4-UC	TZ-600-5-UC
	Transport	Hausted 610	TZ-610-4	TZ-610-4-UC	TZ-610-5-UC
	Transport	Hausted 615	TZ-615-4	TZ-615-4-UC	TZ-615-5-UC
	Transport	Hausted 800	TZ-800-4	TZ-800-4-UC	TZ-800-5-UC
	Transport	Hausted 815	TZ-815-4	TZ-815-4-UC	TZ-815-5-UC
	Horizon Youth Series	Hausted 2000	TZ-2000-4	TZ-2000-4-UC	TZ-2000-5-UC
	Horizon Youth Series	Hausted 4160	TZ-4160-4	TZ-4160-4-UC	TZ-4160-5-UC

NAMC Stretchers BY SIZE

NAMC	Size	Taper	Standard 4"	Ultra Comfort 4"	Ultra Comfort 5"
	24"x76"	8" taper @ head	TZ-76244-D-8	TZ-76244-D-8-UC	TZ-76245-D-8-UC
	26"x76"	8" taper @ head	TZ-76264-D-8	TZ-76264-D-8-UC	TZ-76265-D-8-UC
	30"x76"	8" taper @ head	TZ-76304-D-8	TZ-76304-D-8-UC	TZ-76305-D-8-UC
	24"x76"	12" taper @ head; 4" taper @ foot	76244-B-12-4	TZ-76244-B-12-4-UC	TZ-76245-B-12-4-UC
	26"x76"	12" taper @ head; 4" taper @ foot	76264-B-12-4	TZ-76264-B-12-4-UC	TZ-76265-B-12-4-UC
	30"x76"	12" taper @ head; 4" taper @ foot	76304-B-12-4	TZ-76304-B-12-4-UC	TZ-76305-B-12-4-UC



TZ Medical
Sparked by your ideas

Operating Table Pads






Standard Operating Table Pad

Cost-effective performance and patient comfort

Multi-layered, high resiliency foam for added comfort, support, and durability. Designed for shorter procedures or short stay procedural/exam tables.

2" thick CertiPUR-US foam provides comfortable support for individuals up to 300lbs.

- Soft conforming layer
- Supportive foam base layer



Note: foam and fabric colors may vary

Sealed Seams

- Prevents liquid penetration

Vynmed™ Healthcare Vinyl Cover

- Durable, soft-supple feel offers exceptional comfort
- Anti-bacterial, anti-fungal, water resistant, stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides unique protection; if the outer cover is compromised, foam is protected



Secondary Fluid Barrier



CertiPUR-US Foam



12 Month Warranty



Made in the USA



Vynmed™ Cover





Patient Weight Limit



Multi-Layer Performance Foam

Soft Conforming Layer

Firm Supportive Base



Armboard pads match bed height

Cutout matches the frame opening

Vynmed™ Vinyl Cover



Memory Foam Operating Table Pad

Technology-enhanced for patient comfort

Multi-layered, high resiliency foam for added comfort, support, and durability. Designed for shorter procedures or short stay procedural/exam tables.

3" or 4" thick CertiPUR-US foam provides comfortable support for individuals up to 300lbs.

- Cool Gel infused memory foam top layer
- Performance soft foam internal layer
- Supportive foam base layer

Sealed Seams


- Prevents liquid penetration

Vynmed™ Healthcare Vinyl Cover

- Durable, soft-supple feel offers exceptional comfort
- Anti-bacterial, anti-fungal, water resistant stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides unique protection; if the outer cover is compromised, foam is protected




Note: foam and fabric colors may vary


Cool Gel
Memory Foam

Soft Conforming
Foam

Firm Supportive
Base



- Secondary Fluid Barrier
- CertiPUR-US Foam
- 12 Month Warranty
- Made in the USA
- Vynmed™ Cover
- Patient Weight Limit
- Cool Gel Memory Foam
- Multi-Layer Performance Foam



Armboard pads match bed height

Cutout matches the frame opening

Vynmed™ Vinyl cover

Operating Table Pads

Marathon Mattress	Model	Standard	Ultra Comfort
	Amsco 1080 / 2080 V-Style Mattress (3 piece set)	TZ-5701-2 (2" Thick) TZ-5701-3 (3" Thick)	TZ-5701-3-UC (3" Thick) TZ-5701-4-UC (4" Thick)
	Amsco 3080 V-Style Mattress (3 piece set)	TZ-3080-2 (2" Thick) TZ-3080-3 (3" Thick)	TZ-3080-3-UC (3" Thick) TZ-3080-4-UC (4" Thick)
	Amsco 3085 Mattress (3 piece set)	TZ-3085-2 (2" Thick) TZ-3085-3 (3" Thick)	TZ-3085-3-UC (3" Thick) TZ-3085-4-UC (4" Thick)
	Skytron 5000 / 6000 Mattress (3 piece set)	TZ-6000-2 (2" Thick)	TZ-6000-4-UC (4" Thick)
	Shampaine 4900 Mattress (3 piece set)	TZ-5712-2 (2" Thick)	TZ-5712-4-UC (4" Thick)
	Amsco 3080 Armboard (pair)	TZ-3080-2-AB (2") TZ-3080-3-AB (3") TZ-3080-4-AB (4")	TZ-3080-3-AB-UC (3" Thick) TZ-3080-4-AB-UC (4" Thick)
	Amsco 3085 Armboard (pair)	TZ-3085-2-AB (2") TZ-3085-3-AB (3") TZ-3085-4-AB (4")	TZ-3085-3-AB-UC (3" Thick) TZ-3085-4-AB-UC (4" Thick)
	Skytron 5501 Armboard (pair)	TZ-5501-2-AB (2") TZ-5501-3-AB (3") TZ-5501-4-AB (4")	TZ-5501-3-AB-UC (3" Thick) TZ-5501-4-AB-UC (4" Thick)



X-Ray, Angio & Cath Lab

Designed for those shorter procedures or short stay procedural/exam tables

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Vynmed™ Healthcare Vinyl Cover


- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Phthalate free, latex free and odor resistant

STANDARD X-RAY/ANGIO PADS


- 2" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.
- Soft Conforming foam layer
- Supportive foam base layer

ULTRA COMFORT X-RAY/ANGIO PADS

- 2" or 3" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.
- Cool Gel infused memory foam layer will evenly distribute body weight and cooler sleep experience
- Supportive foam base layer



Note: foam and fabric colors may vary



Note: foam and fabric colors may vary



Sealed Seams



Vynmed™ Vinyl Cover



Patient Weight Limit



Secondary Fluid Barrier




12 Month Warranty



Made in the USA

STANDARD ECHO PADS

- 6" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.
- Soft Conforming foam layer
- Supportive foam base layer
- Includes removable insert to manipulate transducer for optimal visualization



Note: foam and fabric colors may vary

Various shapes available



X-Ray/Angio and Echo Pads

Standard X-Ray/Angio Pads

Model	Order No.
Standard Phillips 134 (95" X 22")	TZ-61-2855
Standard GE Omega (118" X 26")	TZ-61-2745
Standard XRE Angio 2 Table Pad (131.5" X 26.5")	TZ-61-8621
Standard Toshiba Cath Lab Table Pad (85.5" X 28")	TZ-61-2588

Ultra Comfort X-Ray/Angio Pads

Model	Order No.
Secure Core Siemens Koordinat 3 (85" X 29")	TZ-61-2961-UC
Secure Core GE Omega (118" X 26")	TZ-61-2745-UC
Secure Core XRE Angio 2 Table Pad (131.5" X 26.5")	TZ-61-8621-UC
Secure Core Toshiba Cath Lab Table Pad (85.5" X 28")	TZ-61-2588-UC

Standard Echo Pads

Model	Order No.
Standard Echo Pad (75" X 24.5") w/ removable insert	TZ-7524.56-B-12-4-ECHO

Call for additional models



TZ Medical
Sparked by your ideas

Comfort Overlay Pads

Securely attaches over existing mattress






Standard Comfort Overlay Pad

Patients appreciate outstanding comfort
Fits Labor and Delivery and all beds

Designed to attach securely over an existing mattress assuring outstanding comfort and durability with no slipping or shifting.

3" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.

- 3" soft lofted engineered top layer



Nymed™ Healthcare Grade Nylon Top Cover

- Water resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Healthcare Vinyl Bottom Cover

- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area

Easy Storage

- Self storing straps make storing easy and convenient

Birthing Bed shown,
many sizes available.



Patient Weight Limit

Engineered Foam

Attachment Straps

Removable Nylon Cover

Heavy-duty Zipper

Fluid Protection Flap Over Zipper

Secondary Fluid Barrier

CertiPUR-US Foam

12 Month Warranty

Made in the USA




Ultra Comfort Overlay Pad

Patients appreciate outstanding comfort
Fits Labor and Delivery and all beds

Designed to attach securely over an existing mattress assuring outstanding comfort and durability with no slipping or shifting.

3" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.

- Cool Gel infused memory foam top layer
- Performance soft foam internal layer
- Supportive foam base layer



Nymed™ Healthcare Grade Nylon Top Cover

- Water resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Healthcare Vinyl Bottom Cover

- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area

Easy Storage

- Self storing straps make storing easy and convenient

Birthing Bed shown, many sizes available.



Self-storing & Attachment Straps

Nymed 6™ Nylon Top Cover

Vynmed™ Vinyl Bottom Cover

- Cool Gel Memory Foam
- Attachment and Storing Straps
- Removable Nylon Cover
- Patient Weight Limit
- Heavy-duty Zipper
- Fluid Protection Flap Over Zipper
- Secondary Fluid Barrier
- CertiPUR-US Foam
- 12 Month Warranty
- Made in the USA

Comfort Overlay Pads

Patients value exceptional comfort, and our Marathon Overlay delivers just that. Whether it's for extended stays or simply enhancing the comfort of a bed or stretcher, our overlay is the solution. Engineered to be secure, durable, and budget-friendly, it offers a superior experience.

Our Comfort Pads are meticulously crafted to easily attach over existing mattresses, ensuring unparalleled comfort and longevity without any slipping or shifting. The overlay can be conveniently stored using the retention straps, making storage hassle-free. Featuring a top cover made of Nymed™ healthcare nylon and a bottom cover of Vynmed™ healthcare polyurethane coated vinyl, our product provides a resilient, anti-static, stain-resistant, and non-allergenic mattress cover.

Size	Standard Overlay	Ultra Comfort Overlay
81"x33"x3" (Fits Labor and Delivery Model 304)	TZ-26-304	TZ-26-304-UC
76"x36"x3"	TZ-26-2612	TZ-26-2612-UC
76.5"x34.5"x3"	TZ-26-2610	TZ-26-2610-UC
76"x33"x3"	TZ-26-2602	TZ-26-2602-UC
77"x31"x3"	TZ-26-2624	TZ-26-2624-UC
78"x34.5"x3"	TZ-26-2608	TZ-26-2608-UC
78"x36"x3"	TZ-26-2623	TZ-26-2623-UC
79"x36"x3" (Fits Affinity I, II, III, & IV) (Fits Labor and Delivery Models)	TZ-26-2620	TZ-26-2620-UC
80"x35"x3"	TZ-26-80353	TZ-26-80353-UC
80"x36"x3"	TZ-26-80363	TZ-26-80363-UC
76"x26"x3"	TZ-26-76263	TZ-26-76263-UC
76"x30"x3"	TZ-26-76303	TZ-26-76303-UC



TZ Medical
Sparked by your ideas

General Patient and ICU/CCU

Quality to go the distance™



General Hospital Mattress Collection

Step up to the Marathon™ Series Advanced Care Mattresses from our General Hospital Mattress Collection. Engineered with the unique requirements of longer-stay and high-risk patients in focus, our mattresses offer unparalleled support and comfort.

Whether you opt for the Advanced Care™ or the Assure II™, both models are meticulously

designed to address the needs of patients, particularly in preventing pressure ulcers.

Experience the comfort and luxury with our deluxe Advanced Care™ mattress, featuring Vistotech™ gel memory foam for superior support and relaxation. Alternatively, for patients at higher risk, the Assure II™ offers a high loft, convoluted foam top, ensuring optimal comfort and support.

Features	Assure II™	Advanced Care™	Standard Bariatric Assure II (500)	Heavy Duty Bariatric Assure II (1000)	Standard Bariatric Advanced Care (500)	Heavy Duty Bariatric Advanced Care (1000)
Cool Gel Infused Memory Foam		x			x	x
Water Resistant Cover	x	x	x	x	x	x
Secondary Fluid Barrier (protects foam if outer layer is compromised)	x	x	x	x	x	x
Removable Cover	x	x	x	x	x	x
Fluid Protection Flap Over Zipper	x	x	x	x	x	x
Easy to Clean	x	x	x	x	x	x
FirmEdge™ Passive Restraint System	x	x	x	x	x	x
Weight Capacity	300lbs	300lbs	500lbs	1000lbs	500lbs	1000lbs
USA Made Foam	x	x	x	x	x	x
CertiPUR-US Certified Foam	x	x	x	x	x	x



Advanced Care™

General Patient and ICU/CCU Hospital Mattress with Cool Gel Memory Foam

Select the high comfort and unmatched luxury of our deluxe Advanced Care™ pressure reduction mattress for higher risk patients

6" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.

- Cool Gel infused memory foam top layer
- Performance soft foam internal layer
- Supportive firm foam base layer

FirmEdge™ Passive Restraint System

- Provides support on both sides of the mattress during entry and egress

Nymed™ Healthcare Grade Nylon Top Cover

- Water resistant, flame resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Healthcare Vinyl Bottom Cover




- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area



- Cool Gel Memory Foam
- FirmEdge™ Passive Restraint
- Heavy-duty Zipper
- Fluid Protection Flap Over Zipper
- Flame Resistant
- Secondary Fluid Barrier
- CertiPUR-US Foam
- 12 Month Warranty
- Made in the USA
- Removable Nylon Cover
- Patient Weight Limit



Assure II™

General Patient and ICU/CCU Hospital Mattress

Multi-layered, multi-density for optimal comfort and postural support

6" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.

- 2" soft lofted engineered top layer
- Performance soft foam internal layer
- Supportive firm foam base layer

FirmEdge™ Passive Restraint System

- Provides support on both sides of the mattress during entry and egress

Nymed™ Healthcare Grade Nylon Top Cover

- Water resistant, flame resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Healthcare Vinyl Bottom Cover


- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area




Soft Lofted Engineered Foam

Performance Soft Layer

Supportive Base

FirmEdge™ Passive Restraint System



Nymed™ Nylon Top Cover

Fluid Flap Over Heavy-Duty Zipper

Vynmed™ Vinyl Bottom Cover



Heavy-duty Zipper



Fluid Protection Flap Over Zipper



Flame Resistant



Secondary Fluid Barrier



CertiPUR-US Foam



12 Month Warranty



Made in the USA



Removable Nylon Cover



Patient Weight Limit



Multi-Layer Performance Foam



FirmEdge™ Passive Restraint

General Hospital Mattresses

Standard (Weights up to 300lbs) By Size

Size	Assure II	Advanced Care	Size	Assure II	Advanced Care	Size	Assure II	Advanced Care
72"x36"x6"	TZ-42-72366	TZ-44-72366	80"x33"x6"	TZ-42-80336	TZ-44-80336	82"x36"x6"	TZ-42-82366	TZ-44-82366
74"x38"x6"	TZ-42-74386	TZ-44-74386	80"x34"x6"	TZ-42-80346	TZ-44-80346	82"x38"x6"	TZ-42-82386	TZ-44-82386
75"x35"x6"	TZ-42-75356	TZ-44-75356	80"x35"x6"	TZ-42-80356	TZ-44-80356	84"x32"x6"	TZ-42-84326	TZ-44-84326
76"x35"x6"	TZ-42-76356	TZ-44-76356	80"x36"x6"	TZ-42-80366	TZ-44-80366	84"x34"x6"	TZ-42-84346	TZ-44-84346
78"x35"x6"	TZ-42-78356	TZ-44-78356	80"x38"x6"	TZ-42-80386	TZ-44-80386	84"x35"x6"	TZ-42-84356	TZ-44-84356
			82"x35"x6"	TZ-42-82356	TZ-44-82356	84"x36"x6"	TZ-42-84366	TZ-44-84366

Standard (Weights up to 300lbs) By OEM

Model	Assure II	Advanced Care
Stryker Secure II Standard Replacement Mattress	TZ-42-84356	TZ-44-84356
Stryker Acute Care 72" long Standard Replacement Mattress	TZ-42-72366	TZ-44-72366
Stryker Acute Care 80" long Standard Replacement Mattress	TZ-42-80366	TZ-44-80366
Stryker Acute Care 84" long Standard Replacement Mattress	TZ-42-84366	TZ-44-84366
Stryker Go bed 84" long Standard Replacement Mattress	TZ-42-84366	TZ-44-84366
Hill-Rom Advanced 2000 Standard Replacement Mattress	TZ-42-82366	TZ-44-82366
Hill-Rom Advanced 1000 Standard Replacement Mattress	TZ-42-82366	TZ-44-82366
Hill-Rom Centra 850 Standard Replacement Mattress	TZ-42-82366	TZ-44-82366
Hill-Rom model 840 Standard Replacement Mattress	TZ-42-82366	TZ-44-82366
Hill-Rom model 835, 820 and Century Standard Replacement Mattress	TZ-42-80366	TZ-44-80366
Hill-Rom model 820 Standard Replacement Mattress	TZ-42-80366	TZ-44-80366
Hill-Rom Century Standard Replacement Mattress	TZ-42-80366	TZ-44-80366



Bariatric Assure II™

General Patient and ICU/CCU Hospital Mattress

Multi-layered, Multi-density for optimal comfort and postural support. Models supporting weights up to 500lbs or 1000lbs.

7" thick Certi-PUR-US foam provides comfortable support for individual's up to 500 or 1000lbs.

- 2" soft lofted engineered top layer
- Performance soft foam internal layer
- Supportive firm foam base layer

FirmEdge™ Passive Restraint System

- Provides support on both sides of the mattress during entry and egress

Nymed™ Healthcare Grade Nylon Top Cover

- Water resistant, flame resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Healthcare Vinyl Bottom Cover


- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area




Soft Lofted Engineered Foam


Performance Soft Layer

Semi-firm Layer

Supportive Base



FirmEdge™ Passive Restraint System



Nymed™ Nylon Top Cover

Fluid Flap Over Heavy-Duty Zipper

Vynmed™ Vinyl Bottom Cover



Secondary Fluid Barrier



CertiPUR-US Foam



12 Month Warranty



Made in the USA



Removable Nylon Cover



Patient Weight Limit



Multi-Layer Performance Foam



FirmEdge™ Passive Restraint



Heavy-duty zipper



Fluid Protection Flap Over Zipper



Flame Resistant



Bariatric Advanced Care™

General Patient and ICU/CCU Hospital Mattress with Cool Gel Memory Foam. Models Supporting weights up to 500lbs or 1000lbs.

Select the high comfort and unmatched luxury of our deluxe Bariatric Advanced Care™ pressure reduction mattress for your higher risk patients.

7" thick CertiPUR-US foam provides comfortable support for individual's up to 500 or 1000lbs.

- Cool Gel infused memory foam top layer
- Performance soft foam internal layer
- Supportive firm foam base layer

FirmEdge™ Passive Restraint System

- Provides support on both sides of the mattress during entry and egress

Nymed™ Healthcare Grade Nylon Top Cover

- Water resistant, flame resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Vynmed™ Healthcare Vinyl Bottom Cover


- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Latex free and odor resistant

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Removable Cover

- Heavy duty zipper on bottom allows for internal mattress inspection and easy replacement
- Includes fluid protection flap over zipper which directs fluids over and away from the zipper area




Cool Gel Memory Foam

Performance Soft Layer

Semi-firm Layer

Supportive Base

FirmEdge™ Passive Restraint System



Nymed™ Nylon Top Cover

Fluid Flap Over Heavy-Duty Zipper

Vynmed™ Vinyl Bottom Cover



Secondary Fluid Barrier



CertiPUR-US Foam



12 Month Warranty



Made in the USA



Removable Nylon Cover



Patient Weight Limit



Cool Gel Memory Foam



FirmEdge™ Passive Restraint



Heavy-duty zipper



Fluid Protection Flap Over Zipper



Flame Resistant

Bariatric General Hospital Mattresses

Standard Bariatric (Weights up to 500lbs)

Size	Bariatric Assure II	Bariatric Advanced Care
72"x36"x7" 500lbs	TZ-42-HD-72367	TZ-44-HD-72367
80"x34"x7" 500lbs	TZ-42-HD-80347	TZ-44-HD-80347
80"x35"x7" 500lbs	TZ-42-HD-80357	TZ-44-HD-80357
80"x36"x7" 500lbs	TZ-42-HD-80367	TZ-44-HD-80367
80"x39"x7" 500lbs	TZ-42-HD-80397	TZ-44-HD-80397
80"x42"x7" (Wide) 500lbs	TZ-42-HD-80427	TZ-44-HD-80427
80"x48"x7" (Extra Wide) 500lbs	TZ-42-HD-80487	TZ-44-HD-80487
82"x35"x7" 500lbs	TZ-42-HD-82357	TZ-44-HD-82357
84"x34"x7" 500lbs	TZ-42-HD-84347	TZ-44-HD-84347
84"x35"x7" 500lbs	TZ-42-HD-84357	TZ-44-HD-84357
84"x36"x7" 500lbs	TZ-42-HD-84367	TZ-44-HD-84367

Heavy Duty Bariatric (Weights up to 1000lbs)

Size	Bariatric Assure II	Bariatric Advanced Care
72"x36"x7" 1000lbs	TZ-42-HDM-72367	TZ-44-HDM-72367
80"x34"x7" 1000lbs	TZ-42-HDM-80347	TZ-44-HDM-80347
80"x35"x7" 1000lbs	TZ-42-HDM-80357	TZ-44-HDM-80357
80"x36"x7" 1000lbs	TZ-42-HDM-80367	TZ-44-HDM-80367
80"x39"x7" 1000lbs	TZ-42-HDM-80397	TZ-44-HDM-80397
80"x42"x7" (Wide) 1000lbs	TZ-42-HDM-80427	TZ-44-HDM-80427
80"x48"x7" (Extra Wide) 1000lbs	TZ-42-HDM-80487	TZ-44-HDM-80487
82"x35"x7" 1000lbs	TZ-42-HDM-82357	TZ-44-HDM-82357
84"x34"x7" 1000lbs	TZ-42-HDM-84347	TZ-44-HDM-84347
84"x35"x7" 1000lbs	TZ-42-HDM-84357	TZ-44-HDM-84357
84"x36"x7" 1000lbs	TZ-42-HDM-84367	TZ-44-HDM-84367



TZ Medical
Sparked by your ideas



V36™ Deck Bed

Fits the VersaCare® Bed

Compressible foot section adds length, flexibility and accommodates frame expansion/retraction

Dimensions: 86" X 35" X 8"

Retracts to: 75" X 35" X 8"

Model Number: TZ-44-V36

Nymed 6™ Healthcare Grade Stretchable Cast Coated Polyurethane Top Cover:

Designed for patient comfort and clinical efficacy:

- Low friction for ease of patient transfers
- Anti-microbial, anti-static, stain resistant, abrasion resistant
- Non-allergenic, latex free, water resistant

Tri-Zone Foam Construction:

3 independently designed zones:
Head Zone, Seat Zone, Foot Zone

Foam Layers


- Soft plush foam top layer
- Cool Gel infused soft conforming memory foam layer
- Internal supportive middle layer
- Ultra-high density firm base layer

Bottom Cover:

Durable bottom non-skid vinyl cover that is tear resistant and extends mattress performance.

4 built-in carrying handles

that are conveniently located for mattress handling; (includes cleanable coated webbing).



Note: foam and fabric colors may vary



Made in the USA



**All foam
CertiPUR-US
certified**



Side Rail Pads

Models to fit all major brands



Our standard Side Rail Pads attach securely over existing siderails with no shifting and are easily removed when needed.

Vynmed™ Healthcare Vinyl Cover


- Anti-bacterial, anti-fungal, water & stain resistant
- Non-allergenic, non-irritating
- Phthalate free, latex free and odor resistant

CertiPUR-US foam

- Firm foam core offers some protection for patients

Medical Grade Attachment Straps with Stainless Steel Snaps

- Help prevent shifting and ensure the side rail pads are attached securely (depends on models)



Vynmed™ Healthcare Vinyl Cover
Waterproof and tear resistant, polyurethane coated

Firm Foam Core
Soft and supportive

Vynmed™ Healthcare Vinyl Cover
Easily cleaned with mild soap and water

Note: foam and fabric colors may vary



CertiPUR-US Foam



12 Month Warranty



Made in the USA



Vynmed™ Cover

Side Rail Pads for Beds (4 pad models)



Clear vinyl over controls for easy cleaning and protection

Side Rail Pads for Stretchers (2 pad models)



Side Rail Pads

Features	Side Rail Pads for Beds	Side Rail Pads for Stretchers
Sealed Waterproof Cover	x	x
Number of Pads in Set	4	2
Easy to Clean	x	x
Shear Liner	x	x
Latex Free	x	x
Medical Grade Attachment Straps	x	x
CertiPUR-US Certified Foam	x	x
USA Made Foam	x	x

Side Rail Replacement Pads for Beds

Model	Order No.
Stryker GO Bed (4 piece set)	TZ-32-8265
Hill-Rom Advance 2000 (4 piece set)	TZ-2-2000
Hill-Rom Progressa (4 piece set)	TZ-32-P600
Spirit Select (4 piece set)	TZ-32-5700
Stryker Secure II (4 piece set)	TZ-32-8579
Arjo Enterprise 9000 (4 piece set)	TZ-32-E900
Linet Elegana (4 piece set)	TZ-32-7743
Linet Multi-Care (4 piece set)	TZ-32-7761
Stryker in Touch (4 piece set)	TZ-32-8267
Stryker 3002 (4 piece set)	TZ-32-3002

Side Rail Replacement Pads for Stretchers

Model	Order No.
Stretcher Side Rail Pad Full Length (2 each)	TZ-32-7684
Stretcher Side Rail Pad 3/4 Length (2 each)	TZ-32-6084



Labor & Delivery

Stryker Model 304 Beds

Designed for patient comfort and clinical efficacy



Patients will appreciate the exceptional comfort and support of this mattress. You will welcome the durability and performance.

Sealed Seams in high fluid areas.

Nymed 6™ Healthcare Grade Nylon Cover

- Water resistant
- Anti-bacterial, anti-fungal, non-allergenic
- Non-irritating, latex free

Secondary Fluid Barrier

- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Medical Grade Attachment Straps

- Stainless Steel Snaps on foot section to attach pad to bed when required

HEAD SECTION

- CertiPUR-US foam provides comfortable support for individual's up to 300lbs.
- Soft conforming top layer.
- Cool Gel infused memory foam layer will evenly distribute body weight and cooler sleep experience
- Performance foam internal layer
- Supportive firm foam base layer

FOOT SECTION

- Soft conforming top layer
- Supportive firm foam base layer

N NYLON Removable Nylon Cover

300 Patient Weight Limit



Secondary Fluid Barrier

12 Month Warranty


Made in the USA



Nymed 6™ Nylon Cover






Note: foam and fabric colors may vary



Note: foam and fabric colors may vary

Stryker 304 Birthing Bed

Stryker	Size	Type	Order No.
	Stryker 304 Head Section	Head	TZ-304-H
	Stryker 304 Head Section w/Air	Head	TZ-304-H-Air
	Stryker 304 Foot Section	Foot	TZ-304-F



Labor & Delivery

Hill-Rom Affinity™ Series Bed

Designed for patient comfort and clinical efficacy



Your patients will appreciate the exceptional comfort and support of this mattress. You will welcome the durability and performance.

Sealed Seams in high fluid areas.

Nymed 6™ Cover

- Durable, soft-supple feel offers exceptional comfort
- Anti-bacterial, anti-fungal, water resistant, stain resistant
- Non-allergenic, non-irritating
- Phthalate free, latex free and odor resistant

Secondary Fluid Barrier


- Assembled between foam and cover provides a unique protection; if the outer layer is compromised, foam is protected

Medical Grade Attachment Straps

- Stainless Steel Snaps on foot section to attach pad to bed when required

L&D HEAD SECTION


- 5" thick CertiPUR-US foam provides comfortable support for individual's up to 300lbs.
- Cool Gel infused memory foam layer to evenly distribute body weight and cooler sleep experience
- Performance soft foam internal layer
- Supportive foam base layer



Note: foam and fabric colors may vary

L&D FOOT SECTION

- 3.5" thick CertiPUR-US foam, multi-layered
- Soft conforming layer
- Supportive foam base layer



Note: foam and fabric colors may vary

VINYL Nymed 6™ Waterproof Cover

300 Patient Weight Limit

Cool Gel Cool Gel Memory Foam (Head Section)

Attachment Straps







Secondary Fluid Barrier

CertiPUR-US CertiPUR-US Foam















12 Month Warranty

Made in the USA

Hill-Rom Affinity & Affinity II Series

Hill-Rom	Size	Type	Order No.
	Hill-Rom Affinity & Affinity II Straight Head Section	Head	TZ-3613
	Hill-Rom Affinity & Affinity II Straight Head Section w/Air	Head	TZ-3611-AIR
	Hill-Rom Affinity & Affinity II Straight Foot Section - w/folds	Foot	TZ-3612-F
	Hill-Rom Affinity & Affinity II V-Cut Head Section	Head	TZ-3617
	Hill-Rom Affinity & Affinity II V-Cut Head Section w/Air	Head	TZ-3615-AIR
	Hill-Rom Affinity & Affinity II V-Cut Foot Section w/folds	Foot	TZ-3616-F

Hill-Rom Affinity III & IV Series

Hill-Rom	Size	Type	Order No.
	Hill-Rom Affinity III & Affinity IV Straight Head Section	Head	TZ-3627
	Hill-Rom Affinity III & Affinity IV Straight Head Section w/ Air Head	Head	TZ-3627-AIR
	Hill-Rom Affinity III & Affinity IV Straight Double Cover Head Section (includes cover and slip cover)	Head	TZ-3627-DBCVR
	Hill-Rom Affinity III & Affinity IV Straight Head Section w/ Air Head & Double Cover (includes cover and slip cover)	Head	TZ-3627-AIR-DBCVR
	Hill-Rom Affinity III & Affinity IV Straight Head Section - Slip Cover Only	Head	TZ-3627-OVC-CVR
	Hill-Rom Affinity III & Affinity IV Straight Foot Section w/folds & Double Cover (includes cover and slip cover, does not include bottom hook)	Foot	TZ-3626-F
	Hill-Rom Affinity III & Affinity IV Straight Foot Section w/folds - Slip Cover Only	Foot	TZ-3626-OVC-CVR
	Hill-Rom Affinity III & Affinity IV V-Cut Head Section	Head	TZ-3623
	Hill-Rom Affinity III & Affinity IV V-Cut Head Section w/ Air Head	Head	TZ-3623-AIR
	Hill-Rom Affinity III & Affinity IV V-Cut Double Cover Head Section (includes cover and slip cover)	Head	TZ-3623-DBCVR
	Hill-Rom Affinity III & Affinity IV V-Cut Head Section w/ Air Head & Double Cover (includes cover and slip cover)	Head	TZ-3623-AIR-DBCVR
	Hill-Rom Affinity III & Affinity IV V-Cut Head Section - Slip Cover Only	Head	TZ-3623-OVC-CVR
	Hill-Rom Affinity III & Affinity IV V-Cut Foot Section w/folds & Double Cover (includes cover and slip cover, does not include bottom hook)	Foot	TZ-3622-F
	Hill-Rom Affinity III & Affinity IV V-Cut Foot Section w/folds - Slip Cover Only	Foot	TZ-3622-OVC-CVR



Pediatric Mattresses & Pads

Newborn Crib, Infant Pads & Infant Warmer with Support Insert

Sealed Seams in high fluid areas

Vynmed™ Healthcare Vinyl Cover

- Durable, soft-supple feel offers exceptional comfort
- Phthalate Free
- Water and stain resistant
- Anti-bacterial, anti-fungal
- Non-allergenic, non-irritating
- Latex free and odor resistant


1", 1.5" or 2" thick CertiPUR-US foam

- Soft, supportive conforming foam core

Infant Care Stations

- Many infant warmers require support insert

Water/Fluid Proof Newborn Mattress & Pads




Vynmed™ Healthcare Vinyl Cover
Waterproof and tear resistant, polyurethane coated

Conforming Foam
Soft and supportive

Vynmed™ Healthcare Vinyl Cover
Easily cleaned with mild soap and water

Note: foam and fabric colors may vary

Water/Fluid Proof Infant Warmer with Insert




Vynmed™ Healthcare Vinyl Cover
Waterproof and tear resistant, polyurethane coated

Conforming Foam
Soft and supportive

Support Insert

Vynmed™ Healthcare Vinyl Cover
Easily cleaned with mild soap and water

Note: foam and fabric colors may vary



Vynmed™ Waterproof Vinyl Cover



Sealed Seams



Support Insert



Patient Weight Limit Newborn



Patient Weight Limit Infant Warmers



Vynmed™ Vinyl Cover



Waterproof



Easy to Clean



Phthalate Free



CertiPUR-US Foam



12 Month Warranty



Made in the USA



Pediatric Mattresses & Pads

Crib & Youth Mattresses



Sealed Seams

Vynmed™ Healthcare Vinyl Cover

- Durable, soft-supple feel offers exceptional comfort
- Phthalate Free
- Water and stain resistant
- Anti-bacterial, anti-fungal
- Non-allergenic, non-irritating
- Latex free and odor resistant

Crib & Youth Mattress Core

- 3.5" - 6" thick CertiPUR-US foam
- Soft conforming layer
- Supportive foam base layer
- This model does not include plastic support insert



- Sealed Seams
- Support Insert
- Patient Weight Limit
Crib & Youth
- Vynmed™ Vinyl Cover
- Waterproof
- Easy to Clean
- Phthalate Free
- CertiPUR-US Foam
- 12 Month Warranty
- Made in the USA

Vynmed™ Waterproof Vinyl Cover

Breathable Super-A-Vent™

Pediatric Mattresses

Newborn Crib

Size	Order No.
24" x 12" x 1.5"	TZ-27-24121.5-S
24" x 13.5" x 1"	TZ-27-2413.51-S
25" x 12" x 1"	TZ-27-25121-S
25" x 12" x 2"	TZ-27-25122-S
25" x 18" x 1"	TZ-27-25181-S
25" x 20" x 1"	TZ-27-25201-S
25" x 22" x 1"	TZ-27-25221-S
25.5" x 12" x 1"	TZ-27-25.5121-S
25.5" x 12" x 1.5"	TZ-27-25.5121.5-S
26" x 12" x 1.5"	TZ-27-26121.5-S
26" x 12" x 2"	TZ-27-26122-S
26" x 12.5" x 1"	TZ-27-2612.51-S
26" x 12.5" x 2"	TZ-27-2612.52-S
26" x 13" x 2"	TZ-27-26132-S
26" x 13.5" x 2"	TZ-27-2613.52-S
27" x 13" x 2"	TZ-27-27132-S
27" x 16" x 2"	TZ-27-27162-S
28" x 16" x 2"	TZ-27-28162-S
33" x 13.5" x 2"	TZ-27-3313.52-S

Crib & Youth

Size	Order No.	Size	Order No.
44" x 30" x 4"	TZ-39-44304	57" x 27" x 4"	TZ-39-57274
48" x 24" x 4"	TZ-39-48244	57" x 32" x 4"	TZ-39-57324
52" x 27" x 5"	TZ-39-52275	58" x 27" x 4"	TZ-39-58274
52" x 28" x 4"	TZ-39-52284	58" x 32" x 6"	TZ-39-58326
55" x 27" x 5"	TZ-39-55275	58" x 32.75" x 4"	TZ-39-5832.754
55" x 32" x 5"	TZ-39-55325	58" x 34" x 4"	TZ-39-58344
56" x 26" x 3.5"	TZ-39-56263.5	58.5" x 33.5" x 5"	TZ-39-58.533.55
56" x 26" x 3.5" (Pedigo)	TZ-5860002	60" x 26" x 4"	TZ-39-60264
56" x 26" x 5"	TZ-39-56265	67" x 31" x 4"	TZ-39-67314
56" x 27" x 4"	TZ-39-56274	68" x 32" x 6"	TZ-39-68326
56" x 27" x 6"	TZ-39-56276	69" x 32" x 4"	TZ-39-69324
56" x 31" x 4"	TZ-39-56314	69" x 33" x 4"	TZ-39-69334
56.5" x 26.5" x 4"	TZ-39-56.526.54	70" x 33" x 5"	TZ-39-70335

Many other sizes available
please call for part number

Pediatric Mattress Pads

Infant Pads

Size	Order No.
22" x 18" x 1"	TZ-27-22181-S
22" x 18" x 1" (pack of 5)	TZ-27-22181
23" x 14.5" x 1"	TZ-27-2314.51-S
23.5" x 17.5" x 1"	TZ-27-23.517.51-S
24" x 13" x 1"	TZ-27-24131-S
24" x 13" x 1" (pack of 5)	TZ-27-24131
24" x 18" x 1"	TZ-27-24181-S
24" x 18" x 1.5"	TZ-27-24181.5-S
24.5" x 13.5" x 2"	TZ-27-24.513.52-S
25" x 13.5" x 2"	TZ-27-2513.52-S
25" x 14" x 1"	TZ-27-25141-S
25" x 17" x 1"	TZ-27-25171-S
25" x 18" x 1"	TZ-27-25181-S
25" x 18" x 1" (pack of 5)	TZ-27-25181
25" x 18.5" x 1.5"	TZ-27-2518.51.5-S
25.5" x 12" x 1"	TZ-27-25.5121-S
25.5" x 13.5" x 1"	TZ-27-25.513.51-S
25.5" x 13.5" x 1" (pack of 5)	TZ-27-25.513.51
25.5" x 13.5" x 2"	TZ-27-25.513.52-S
25.5" x 13.5" x 2" (pack of 5)	TZ-27-25.513.52
26" x 19" x 1"	TZ-27-26191-S
26" x 20" x 1"	TZ-27-26201-S
26" x 20" x 1" (pack of 5)	TZ-27-26201
26" x 20.5" x 1"	TZ-27-2620.51-S
26" x 22" x 1"	TZ-27-26221-S
27.5" x 14" x 2"	TZ-27-27.5142-S
28" x 17.5" x 1" (Chaparral)	TZ-27-2817.51-S
28.5" x 22.5" x 1"	TZ-27-28.522.51-S
28.5" x 22.5" x 1.5"	TZ-27-28.522.51.5-S
32" x 23.75" x 1"	TZ-27-3223.751-S
36" x 19" x 1"	TZ-27-36191-S
36" x 19" x 1" (pack of 5)	TZ-27-36191

Infant Warmer with Support Insert

Size	Order No.
23" x 20" x 1"	TZ-28-23201
24" x 14" x 1"	TZ-28-24141
24.5" x 13.5" x 1"	TZ-28-24.513.51
25" x 18" x 1"	TZ-28-25181
25" x 20" x 1"	TZ-28-25201
25.75" x 20.5" x 1"	TZ-28-25.7520.51
26" x 19" x 1"	TZ-28-26191
26" x 20" x 1"	TZ-28-26201
26" x 21" x 1"	TZ-28-26211
30" x 23" x 1"	TZ-28-30231
36" x 24" x 1"	TZ-28-36141

Part numbers are single pads unless noted as pack of 5