

2021-2022 Product Catalog



SHARE YOUR IDEA



PARTNER WITH TZ



EARN ROYALTIES

INNOVATION STARTS WITH YOU

www.tzmedical.com



TZ Medical
Sparked by your ideas



TZ Medical

Sparked by your ideas

TABLE OF CONTENTS

The TZ Medical Story	_____	3
Sparked By Your Ideas	_____	4
Radial Access	_____	5
Femoral Access	_____	12
Patient Positioning	_____	17
Interventional	_____	29
Defibrillation Electrodes	_____	36
Safety & Efficiency	_____	40
Electrophysiology	_____	43
Radiation Protection	_____	47
Procedure Lighting	_____	52
Cardiac Monitors	_____	56
ClearBeat Monitoring	_____	59
Item Number Index	_____	61

TZ Medical Inc.

17750 SW Upper Boones Ferry Rd. #150
Portland, OR 97224
Office: 800.944.0187 Fax: 503.639.0239
www.tzmedical.com

THE TZ MEDICAL STORY

When TZ Medical founder, Tom Tribou, sat planning his new company in a Portland coffee shop in 1990, he made two bedrock decisions. The first was that TZ Medical would turn the paradigm of medical equipment manufacturing on its head by taking the time to listen to medical professionals on the front lines, and find out what they needed, before huddling with its own engineers. The second decision was that pricing should be at the value level, to ensure that these bright ideas would make their way into the hands of critical-care medical professionals faster.

Over the years these ideas have served the company well. What was once a small collection of critical-care devices has grown into a whole family of medical products, sparked by your ideas. With products sold across the nation and all over the world, TZ Medical has built a solid reputation for itself on the first two decisions laid out in the company's infancy.

Today, the TZ Medical team is recognized as a group that creates products based on input from professionals in some of the finest medical facilities throughout the world. These innovative health professionals have often helped supply the missing piece of a puzzle, revealed the essence of a problem, or posed the perfect "what if" question in order to help solve a dilemma, and facilitate practical solutions.

Between our dedication to listening, and our firm belief in affordable pricing, it's no wonder that TZ Medical has been, and will continue to be, a lasting leader in the medical device manufacturing world.



OUR MISSION STATEMENT

Our mission is to design, manufacture, and deliver products that meet and exceed the needs and expectations of each customer. We will sustain our vision and mission by constantly seeking renewal via continuous education and learning, and the application of new technologies and best business practices. We will do business as model individual and corporate citizens in our local and international communities. We will "play to win", acting with absolute honesty, integrity, and fairness in the way we conduct our business and the way we live our lives.



INNOVATION STARTS WITH YOU



SHARE YOUR IDEA



PARTNER WITH TZ



EARN ROYALTIES

SPARKED BY YOUR IDEAS

At TZ Medical, we take pride in our practice of partnering with physicians, nurses, and technologists to help bring new product ideas to market. With decades of practical experience in design, manufacturing, and engineering, we excel in transforming your innovative ideas into product realities. Our most reliable and valuable asset to improving healthcare is YOU.

FIND OUT MORE



800.944.0187



tzmedical.com



TZ Medical
Sparked by your ideas



TZ Medical
Sparked by your ideas

RADIAL ACCESS


With a variety of innovative solutions to choose from, TZ Medical's radial access line of products help optimize patient positioning, hemostasis, and ergonomics to improve clinical outcomes and patient and staff safety.

ARC ● Cobra Board ● StandTall™ ● Flat Radial Board
Secure Radial Board ● Radial Runway

ADJUSTABLE RADIAL CUFF



TZ Medical
Sparked by your ideas




Item Number	Band Size	Diameter	UOM
-------------	-----------	----------	-----

RCB-0047	Small	5-6.5 inches	Box of 5
RCB-0048	Medium	6-7.25 inches	Box of 5
RCB-0049	Large	6.75-8.5 inches	Box of 5



COBRA BOARD

The Cobra Board is a buttress device that can prevent arm drift during left radial artery procedures. This device accommodates a wide range of anatomies due to the flexible 12" gooseneck design. The gooseneck is also high tensile to prevent the instrumented arm from drifting away from the operator during catheterization procedures. The carbon fiber base plate fits between the procedural table and mattress, and the cuff is positioned behind the patient's tricep, for optimal patient comfort and enhanced operating ergonomics.



Features

- High tensile 12" gooseneck design
- Carbon fiber base plate
- Comfortable and ergonomic cuff
- Prevents arm movement during procedures



Item Number	Description	UOM
CBR-0004	Adjustable left arm positioner	1 Each

STANDTALL™



TZ Medical
Sparked by your ideas

DESIGNED TO OPTIMIZE LEFT RADIAL ACCESS

Take advantage of the benefits of Left Radial Artery Access from a comfortable working position. The StandTall™ is an intuitively designed external vascular sheath extender that conveniently redirects workflow to the right side of the procedure table. As a result, you can operate with your preferred dominant hand, maintain optimal posture, and reduce musculoskeletal stress and fatigue.



Physician benefits

- Decreased fatigue and musculoskeletal stress
- Improved site access and stability
- Streamlined workflow
- Allows for a more ergonomic posture for the physician from the right side of the table

Patient benefits


- Greater comfort with left arm positioning
- Quicker recovery and reperfusion times
- Enables convenient access to non-dominant hand (90% of patients are right handed)
- Improved overall procedure satisfaction

Hospital benefits

- Reduced cost vs. femoral approach
- Reduced post-op patient management
- Fewer physician work-related absences due to MSK stress or injury



Without StandTall™



With StandTall™



Item Number	Description	UOM
RD-ST0010L	External vascular sheath extender 10 cm length	Box of 5
RD-ST0015L	External vascular sheath extender 15 cm length	Box of 5
RD-ST0025L	External vascular sheath extender 25 cm length	Box of 5



FLAT RADIAL BOARD

The Flat Radial Board is a lightweight carbon fiber board specifically shaped to allow convenient radial access in conjunction with medical imaging and patient comfort. By combining radial access with a landing zone for equipment, the Flat Radial Board takes radial access to a whole new level.

The Flat Radial Board can be used for radial access on either the right or left side.



The dimensions of the Flat Radial Board are 44.4" x 22.9", offering plenty of space and positioning options.

Features

- Can be used on either the right or left side for radial access
- Large landing area for equipment
- Lightweight carbon fiber construction (No foam core)
- Invisible under fluoroscopy
- Can be used for all diagnostic and interventional imaging procedures

The Flat Radial Board slides under the mattress for increased stability.



Item Number

Description

UOM

CZ-400

Flat Radial Board

1 Each

SECURE RADIAL BOARD



TZ Medical
Sparked by your ideas

The TZ Medical Secure Radial Board is designed to secure and orientate the patient's right or left arm during any vascular procedures. The arm board and straps eliminate unwanted movement yielding a consistent image for the physician when performing vascular access.



Disposable Straps
Two 2" x 24"
One 2" x 18"

Features

- Can be used on either the right or left side for radial access
- Straps reduce patient movement
- Large landing area for equipment
- Lightweight, ultra-durable carbon fiber construction
- Invisible under fluoroscopy
- Can be used for all diagnostic and interventional imaging procedures



Item Number	Description	UOM
CZ-400-TVA	Secure Radial Board with cut outs for straps	Each
TVA-MC-2	Disposable Straps - two 24" straps, one 18" strap	Each (3 straps)




The Radial Runway provides a radiolucent platform that keeps the patient in normal anatomical position for ideal radial access. Disposable straps allow the wrist to be positioned in multiple neutral or extended angles optimizing patient comfort, and the foam pad may be removed from the base at the conclusion of the case. The disposable brace can then be attached to the foam pad, providing support for the patient's wrist during recovery.



Features

- Can be used for radial access on either wrist
- Adjustable straps fit any sized patient
- Disposable straps, foam pad, & brace
- Base can be sterilized and reused
- Keeps patient comfortable during procedure
- Detachable pad with brace remains with patient for recovery
- Anti-skid surface



Item Number	Description	UOM
RRU-0008	Radial Runway - Reusable Base (Non-Sterile)	1 Each
RRU-0009	Radial Runway - Disposable Top Pad (Non-Sterile)	Box of 10



TZ Medical
Sparked by your ideas

FEMORAL ACCESS


TZ Medical's line of femoral access products offers a variety of configurations to accommodate multiple patient and procedure types. Our hemostasis products provide our customers and their patients with time-tested efficiency and unparalleled safety.

Neptune • EZ Hold • Femoral Pressure Discs




The TZ Medical Neptune Pad is a natural product derived from sustainably harvested seaweed that utilizes calcium ions to facilitate hemostasis. These ions are activated by the introduction of blood to the Neptune pad and help to speed up clot formation at the puncture site naturally.

Neptune has no known contraindications and is safe to use on all patients. The texture of the pad eliminates sliding, sliming, and sticking to gloves, and because the entire pad is the active ingredient, there is no "right-side-up."



Features

- Speeds up Hemostasis
- Extra-large surface area
- Naturally derived from seaweed
- Ad-Skin Patch Transparent Film Dressing included
- Latex-free & contains no animal products
- Texture helps prevent sliding or sticking
- Cost-effective




Neptune's calcium alginate provides additional calcium (Ca^{++}) ions for the clotting process and to speed topical hemostasis.





STANDARD NEPTUNE

At 2"x 2", Neptune is one of the largest hydrophilic pads available on the market. Calcium ions introduce elements necessary to the clotting cascade, which helps speed up hemostasis and form a stable clot at the puncture site.



NEPTUNE PLUS

The same great hemostasis results from our standard Neptune Pad are also achieved with Neptune Plus, with the additional benefits of a protective polymer backing. The calcium alginate provides the calcium ions necessary for the hemostasis process and helps to speed up clot formation at the puncture site naturally. The Neptune Plus backing is translucent, allowing for visualization of hemostasis.



SILVERSTAT

The SilverStat hemostasis pad is an all-natural calcium alginate pad infused with silver nitrate. These compounds combine to make a safe and highly effective hemostatic dressing. Silver is known to be one of the most effective agents in preventing bacterial infections in a wound. The element naturally eliminates bacterium, even in low doses.

Item Number	Description	UOM
8870-01	Neptune 2x2 hemostatic pad	Box of 10
8870-03	Neptune Plus 2x2 hemostatic pad w/ polymer backing	Box of 10
8870-07	SilverStat 2x2 antimicrobial hemostatic pad w/ silver nitrate	Box of 10




EZ HOLD & DISCS

The EZ Hold is a manual-assist compression device developed to aide the clinician in manual femoral pressure.

The device was designed and patented by an RN who sought to reduce hand fatigue and direct contact with bodily fluids. The ergonomic design of the EZ Hold permits the clinician to hold pressure comfortably with one hand while palpating the pulse with the other.

The EZ Hold is best used in conjunction with our line of Neptune Discs. These come in two styles, and have the advantage of having a Neptune pad built into the disc, which speeds up hemostasis and reduces holding times.




Features

- Reduces hand fatigue
- Ergonomic design
- Cuts down on direct contact with bodily fluids
- Can be held comfortably with one hand
- Fits a wide variety of discs
- Neptune Discs speed up hemostasis

Item Number	Description	UOM
-------------	-------------	-----

8260-20	2 lb. EZ Hold	1 Each
---------	---------------	--------



COMPRESSION DISCS




TZ Medical
Sparked by your ideas

TZ Medical's Femoral Pressure Discs are used as a hemostasis device with either our EZ Hold or compression clamp products. Our discs help to apply firm and even pressure on the access point to ensure hemostasis is achieved safely and effectively. Our discs come in a variety of shapes and sizes to accommodate all technician preferences and patient body types.

REGULAR DISCS



Round Universal



Oval Universal



Round Omega



Large Round Omega



Oval Omega

NEPTUNE DISCS



Round Universal Neptune



Oval Universal Neptune

Item Number	Description	UOM
8880-01	Oval Universal disc	Box of 50
8880-02	Round Universal disc	Box of 50
8880-03	Oval Universal disc w/ Omega opening	Box of 50
8880-04	Round Universal disc w/ Omega opening	Box of 50
8880-06	Large Round Universal disc w/ Omega opening	Box of 50
8880-07	Oval Universal disc w/ Neptune	Box of 10
8880-08	Round Universal disc w/ Neptune	Box of 10



TZ Medical
Sparked by your ideas

PATIENT POSITIONING

TZ Medical's line of patient positioning products offers a variety of options to safely and effectively position patients for any type of hospital procedure. Whether you need replacement table pads, positioning aids, immobilizers, or need to secure excess patient tissue, we have you covered.

Patient Immobilizer ● PRS ● Comfort Halo
Sidewinder Board ● Comfort Cradle ● SoftGels



The Patient Immobilizer allows the patient to be secured during imaging procedures in order to limit excess movement, thereby eliminating motion artifact. This greatly reduces the need for repeat subtraction runs, which in turn decreases patient and staff radiation exposure.

The patient's wrists can be safely attached to the leg wrap by using the quick snap and rapid adjustment webbing, which can also be quickly undone as soon as the procedure is complete. For additional safety, a second wrap may also be secured around the patient's chest.



Features

- Improves safety for patient and staff during imaging procedures
- Eliminates motion artifact caused by patient movement
- Cuts down on the need for repeat subtraction runs.
- Quick snap and rapid adjustment webbing for fast immobilization of patient
- Patient may be effectively secured during procedures where conscious sedation is used

The Patient Immobilizer XXL Extender allows patients of any size to be safely secured during imaging procedures. The XXL Extender attaches easily to the Patient Immobilizer and offers a much greater range of size adjustment.



Item Number	Description	UOM
CZ-106	Patient Immobilizer (includes 1 set of CZ-008)	1 Kit
CZ-107	GE Patient Immobilizer (includes 1 set of CZ-008)	1 Kit
CZ-008	Disposable Wrist Wrap	Box of 10
CZ-012	XXL Extender	1 Each



2 Pad, 1 Strap System

The TZ Medical Pannus Retention System (PRS) is a hook & loop pad and strap system designed to safely and efficiently retain excess tissue for high BMI patients. It can be used as a patient positioning tool for a variety of procedural tasks, including gaining femoral access for heart caths and TAVR's, deflecting the left arm to access the left radial, and deflecting tissue to access pacemaker/ICD implant pockets.



2 - 5.5" x 7.5" Pads



1 - 96" Strap

Intended for skin applications, the proprietary adhesive in the pads is residue-free and minimizes the risk of skin tears and degradation typically seen with the use of tape and competing products. Additionally, the product can be adjusted after application, giving clinicians complete control over how their patients are positioned.

Features

- Hypoallergenic, latex & residue free
- Adhesive pads will not tear or damage skin
- Ideal for groin access, left side radial access, and more!
- Completely adjustable system - will fit any sized patient
- Can be applied in less than 1 minute
- Easy to reposition after application
- No additional products are required



Item Number	Description	UOM
CZ-PRS-02	Pannus Retention System 2 pad, 1 strap system	Box of 10




4 Pad, 2 Strap System

The TZ Medical Pannus Retention System is used to safely and effectively position and support excess tissue, allowing easy abdominal access during cesarean sections. Easily applied in one minute, the Pannus Retention System is designed to maintain patient comfort while allowing patient positioning from any angle. The device consists of four adhesive pads and two medical-grade hook & loop straps that can be adjusted after application. Intended for skin applications, the proprietary adhesive in the pads is residue-free and minimizes the risk of skin tears and degradation typically seen with the use of tape and competing products.



4 - 5.5" x 7.5" Pads



2 - 48" Strap



Safely retains large amounts of tissue



Features

- Hypoallergenic, latex & residue free
- Adhesive pads will not tear or damage skin
- Ideal for High-BMI Cesarean Sections and Trans Aortic Valve Replacements (TAVR)
- Completely adjustable system - will fit any sized patient
- Can be applied in less than one minute
- Easy to reposition after application
- Sterile and designed for surgical procedures

Item Number

Description

UOM

CZ-PRS-04

Pannus Retention System
4 pad, 2 strap system, sterile

Box of 10



COMFORT HALO

The TZ Medical Comfort Halo is designed to enhance safety and comfort for patients undergoing conscious sedation.

The polycarbonate base is compatible with all bed types, allowing it to easily slide underneath the imaging table mattress pad. The steel and vinyl tubing can then be manipulated into a multitude of different draping positions, thereby keeping the drape out of the patient's face and allowing complete visualization of the patient at all times.



For All Procedures Using Conscious Sedation Including:

- PPM/ICD/Bi-V ICD implants
- Chemo/Port placements
- CVP and PICC line insertion
- Antegrade femoral angiography
- Rt Heart Caths from Subclavian Approach

Features

- Increases patient comfort by keeping anesthesia drape off of patient's face during conscious sedation
- Ideal for all departments using radiographic imaging and sterile technique
- Offers an unobstructed view of patient during procedure
- The flexible design may be molded to the desired shape or height
- Visual reminder to perform "Time Out", enhancing patient safety/ improving staff compliance with JCAHO requirements
- Compatible with any bed type, no arm rails necessary



Item Number	Description	UOM
-------------	-------------	-----

CZ-009

Comfort Halo -
Anesthesia screen

1 Each

SIDEWINDER BOARD



TZ Medical
Sparked by your ideas

The TZ Medical Sidewinder Board is a simple solution that improves patient safety and comfort by effectively positioning the anesthesia drape away from the patient's face while undergoing conscious sedation. Designed with a sturdy carbon fiber base and a flexible gooseneck, the Sidewinder Board slides underneath the table mattress on either side of the patient and is easily adjustable during any procedure. The malleable steel and vinyl gooseneck can be molded into unlimited draping positions, enhancing visualization and patient satisfaction.



Features


- Improves patient comfort
- Enhances visualization of the patient during procedures
- Flexible gooseneck can be molded to any desired shape or height
- Carbon fiber base can be positioned on either side of the patient
- Compatible with any bed type, no arm rail necessary

Item Number	Description	UOM
-------------	-------------	-----

CZ-009-HH	Anesthesia drape positioner	1 Each
-----------	-----------------------------	--------



COMFORT CRADLES



**Polypropylene
Comfort Cradle**

**Original
Comfort
Cradle**



The TZ Medical Comfort Cradle is a versatile product that aids both the patient and the operator during procedures. The cradle's extended width improves comfort, making long procedures more tolerable for the patient. It is available in a clear polycarbonate material or a white polypropylene material, both providing durability and comfort.

The dimensions of the Comfort Cradle provide superior support during vascular procedures when compared with smaller, less robust arm boards.

The Comfort Cradle also offers the ability to position the Flat Detector (or II) in direct contact with the patient in all views, reducing secondary and scatter radiation. The result is superior image quality and reduced radiation exposure for the staff.

Features

- Patient is more accessible to the physician
- Extra wide dimensions for added comfort
- Reduces secondary and scatter radiation
- Options of either durable polycarbonate or polypropylene construction
- Superior image quality
- Comes with a six month warranty



23

Item Number	Description	UOM
CZ-100	Comfort Cradle Kit	Set of 2
CZ-001	Comfort Cradle Replacement Board	1 Each
CZ-100-PP	Polypropylene Comfort Cradle Kit	Set of 2
CZ-001-PP	Polypropylene Comfort Cradle Replacement Board	1 Each
CZ-002	Comfort Cloud Arm Board Padding	Box of 20

SOFTGELS



TZ Medical
Sparked by your ideas

The TZ Medical SoftGels provide superior cushioning and positioning, ensuring optimal exposure to the surgical site.



Features

- Reduces skin changes and pressure induced soreness
- Allows for correction to airway management
- Increases patient positioning stability
- Washable and reusable for any surgical procedure
- Allows for easy exposure to surgical site, providing optimal physiological safety for the patient





TZ-GELDA



TZ-GELDA-3



TZ-GELPD



TZ-GELND



TZ-GELHRS



TZ-GELPHR



TZ-GELSHP



TZ-GELHDP

Item Number	Description
TZ-GELDA	Adult Gel Donut 8"OD x 3"ID x 1.75"H
TZ-GELDA-3	Extra Deep Donut 8"OD x 3"ID x 3"H
TZ-GELPD	Pediatric Donut 5.5"OD x 2.25"ID x 1.25"H
TZ-GELND	Neonatal Donut 3.07"OD x 1.68"ID x 0.81"H
TZ-GELHRS	Adult Horseshoe Donut 8"OD x 3"ID x 1.75"H
TZ-GELPHR	Pediatric Horseshoe Donut 5.5"OD x 2.25"ID x 1.25"H
TZ-GELSHP	Supine Headrest 11"OD x 9.5"ID x 2.625"H
TZ-GELHDP	Head Pad 11"OD x 9.5"ID x 3.5"H





Item Number	Description
TZ-GELAF	Flat Armboard Gel 24" x 4.5" x 0.625"
TZ-GELAC	Gel Contoured Armboard 20" x 5.75" x 2.25"
TZ-GELUPL	Ulnar Nerve Protector, Large 20.13" x 13.75" x 0.25"
TZ-GELUPS	Ulnar Nerve Protector, Small 16.25" x 6.5" x .5"



Item Number	Description
TZ-GELCC	Gel Calf Cup
TZ-GELHC	Gel Heel Cup
TZ-GELHP	Heel Protector Pad
TZ-GELAS	Ankle Stirrup Pad



TZ-GELCC



TZ-GELHC



TZ-GELHP



TZ-GELAS



TZ-GELPP10



TZ-GELPP14



TZ-GELPP18



TZ-GELCR4-12



TZ-GELCR4-16



TZ-GELCR4-20



TZ-GELCR6-12



TZ-GELCR6-16



TZ-GELCR6-20

Item Number	Description
TZ-GELPP10	Gel Patient Positioner 10" x 4.5" x 2" (Sandbag Style)
TZ-GELPP14	Gel Patient Positioner 14" x 4.5" x 2" (Sandbag Style)
TZ-GELPP18	Gel Patient Positioner 18" x 4.5" x 2" (Sandbag Style)
TZ-GELCR4-12	Flat Bottom Gel Chest Roll 12" x 4" x 4"
TZ-GELCR4-16	Flat Bottom Gel Chest Roll 16" x 4" x 4"
TZ-GELCR4-20	Flat Bottom Gel Chest Roll 20" x 4" x 4"
TZ-GELCR6-12	Flat Bottom Gel Chest Roll 12" x 6" x 6"
TZ-GELCR6-16	Flat Bottom Gel Chest Roll 16" x 6" x 6"
TZ-GELCR6-20	Flat Bottom Gel Chest Roll 20" x 6" x 6"





Item Number	Description
TZ-GELFT	Full Length Table Gel Overlay 20" x 72" x 0.5"
TZ-GELFT3-FT	Gel Foot Table Overlay 20" x 20" x 0.5"
TZ-GELFT3-HD	Gel Head Table Overlay 20" x 10" x 0.5"
TZ-GELFT3-T	Gel Torso Table Overlay 20" x 42" x 0.5"
TZ-GELFT3	3 Piece Table Gel Overlay (Includes: TZ-GELFT3-FT, TZGELFT3-HD, TZ-GELFT3-T)



TZ-GELFT



TZ-GELFT3



TZ-GELFT3-FT





TZ Medical
Sparked by your ideas

INTERVENTIONAL

TZ Medical's line of interventional products combines innovation with efficiency to deliver high-quality solutions for interventional procedures. With multiple options to accommodate a variety of procedure types, our products are designed to optimize ergonomics, visualization, and accuracy to improve clinical outcomes for the patient and staff.

AccuFlate ● M-Cath Microcatheter ● Control
Mechanical Thrombectomy System ● Shoot 'N See

ACCUFLATE INFLATION DEVICE



TZ Medical
Sparked by your ideas

The TZ Medical AccuFlate is a 20ml disposable device designed to accurately control and monitor pressure during inflation and deflation of dilation balloon catheters. Our patented technology allows the Inflation Device to generate pressure up to 30 atm with accuracy of ± 1 atm.



Patented self-locking mechanism for easy one-handed locking control, as well as gradual and rapid pressure release ergonomic handle.



Precision manometer with a maximum pressure of 30 atm.



Optional, male rotating luerlock, 3-way stopcock with "ON" or "OFF" tap handle.

Item Number	Description	UOM
-------------	-------------	-----

IDD-0001

Balloon Inflation Device


Box of 5




M-CATH™

Microcatheter

- Small entry profile: 1.6 Fr
- Unique shaft design
- High push proximal shaft
- Enhanced torquability
- High trackability



Tapered and robust tip for enhanced crossability



Hybrid shaft design for maximum pushability and torquability




Ergonomic design for high comfort and control

MICROCATHETER



TZ Medical
Sparked by your ideas



Technical Information

Useable Length	Entry Profile	Dist. Outer Diameter	Max. Outer Diameter	Dist. Inner Diameter	Prox. Inner Diameter	Max. Guidewire Size
135 cm	1.6 Fr'	2.25 Fr	3.3 Fr	0.016"	0.018"	0.014"

'Entry profile average value

Ordering Information

Product Number	Description	UOM
M-CATH	Microcatheter	1 each

Manufactured by: Acrostak (Schweiz) AG, Stegackerstrasse 14, CH-8409 Winterthur, Switzerland.
Contact: sales@acrostak.com T: + 41 (0)22 347 68 70 F: + 41 (0)52 233 95 50 www.acrostak.com
M-CATH is trademark of Acrostak (Schweiz) AG

Caution: Federal law restricts this device to sale by or on the order of a physician



The Control Mechanical Thrombectomy Aspirator combines simplicity with maximum pressure to provide optimal control during thrombectomy procedures. Designed with dual valve technology, vacuum pressure is easily controlled by squeezing and releasing the handle, allowing for continuous force while preventing re-injection of the aspirant. The closed fluid system drains the aspirant directly into the 275 ml collection bag when releasing the handle, providing successful thrombectomy with a fraction of the blood loss compared to pump based systems. Instantly start, stop, increase, decrease, and maintain pulse force with a simple to use system that is easily accessible on the procedure table.



MORE *Power Control*

LESS *Complexity Blood Loss*

Features


- Instantly start, stop, increase, or decrease force
- Complete control of consistent pressure for a fraction of the blood loss
- Reset force with each squeeze-release cycle
- More vacuum power compared to common pump-based systems
- Maintain continuous aspiration force
- 275 ml collection bag
- Dual valve technology creates a completely closed fluid system
- Compatible with all thrombectomy catheters on the market





MECHANICAL THROMBECTOMY SYSTEM



TZ Medical
Sparked by your ideas



HIGH-PERFORMANCE 5F - 11F OTW & RAPID-EXCHANGE THROMBECTOMY CATHETERS





Item Number	Description	Crossing/Shaft	Sheath Size	Guide Wire	Working Length
ASP-1001	Control RX Mechanical Thrombectomy	0.020" / 0.054"	6F	0.014	136 cm
ASP-1005	Control 5F Mechanical Thrombectomy	0.065" / 0.065"	5F	0.035	135 cm
ASP-1009	Control 6F Mechanical Thrombectomy	0.075" / 0.075"	6F	0.035	135 cm
ASP-1013	Control 7F Mechanical Thrombectomy	0.050" / 0.090"	7F	0.035	90 cm
ASP-1021	Control 11F Mechanical Thrombectomy	0.050" / 0.137"	11F	0.035	90 cm
ASP-1029	Control Mechanical Aspirator	Compatible with any thrombectomy catheter			



SHOOT 'N SEE

The TZ Medical Shoot 'N See is a flexible medical grade tape for use during any vascular procedure involving radiography or fluoroscopy. Applied to the skin or drape before the procedure, the radiopaque markings are visible both to the eye and under fluoroscopy, which allows for quick cross-referencing between the inside and outside of a patient's body.

The flexible tape follows patient contours and does not tear under normal use. The Shoot 'N See reduces the need for exposure to contrast media. It also helps with more accurate lesion locating, and more precise stent and catheter sizing.



Features

- Fully visible with or without fluoroscopy
- Arrow indicators to show precise location of lesion(s)
- Allows for quick cross-referencing between inside/outside of patient's body
- Packaged sterile and measures up to 40 cm
- Great for vascular procedures

Item Number	Description	UOM
-------------	-------------	-----

RCM-0005-10	Shoot 'N See 40 cm tape	Box of 10
RCM-0005-30	Shoot 'N See 40 cm tape	Box of 30
RCM-0005-50	Shoot 'N See 40 cm tape	Box of 50



DEFIBRILLATION ELECTRODES

Our multi-function electrodes are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Our full line of electrodes direct-connect to all major brands of defibrillators and include radiotransparent, radiotranslucent, sterile and pediatric options.

Multifunction Adult Electrodes • Pediatric Electrodes
Sterile Electrodes



TZ Medical Defibrillation Electrodes are constructed to adhere and conform to the patient's body, providing superior current distribution, which results in better clinical outcomes.

TZ Medical's Pediatric Defibrillation Electrodes boast the same clinical benefits as our adult pads, however, are designed for your adolescent patients. Our Pediatric Defibrillation Electrodes are constructed to adhere and conform to the patient's body, providing superior current distribution, which results in better clinical outcomes.

Our multi-function electrodes are designed for EP and Cath Labs, Cardiac Surgery, Critical-Care Areas and other related departments. Our multi-function sterile electrodes are ideal for a variety of procedures, including TAVR, AVR, Mini-Maze, ICD Implants, Laser Lead Extraction and much more. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Compatible with biphasic and monophasic defibrillators, our electrodes offer direct connections to a variety of defibrillators, such as Zoll®, Medtronic® QUIK-COMBO™, and Philips®. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.



Features


- Lead wire lengths available in 42", 60" and 120" for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Unique edge effect technology disperses energy uniformly, reducing burns
- Radiotransparent and radiotranslucent options available



Pad Placement Guide



TZ Medical
Sparked by your ideas



Lead Wire Length


TZ Sterile P-121-T1 & P-122-T1: 120"

TZ Electrodes: 60"

TZ Sterile P-224-T1: 42"

Competitor Electrodes: 42"

Competitor Electrodes: 24"





**Medtronic® M1
Connector**



**Philips® P2
Connector**



**Universal T1
Connector**



**Zoll® Z1
Connector**



Product Description	Size	Medtronic® ¹	Philips® ²	Universal ³	Zoll® ⁴
Adult Radiotranslucent pads, leads out	A	P-111-M1	P-111-P2	P-111-T1	P-111-Z1
Adult Radiotransparent pads and leadwires	A	P-211-M1	P-211-P2	P-211-T1	P-211-Z1
Adult Radiotransparent pads and leadwires, leads out	A	P-311-M1			P-311-Z1
Adult Radiotranslucent, sterile	A	P-121-T1 PAD 1111 Adapter*	P-121-T1 PAD 109 Adapter*	P-121-T1	P-121-T1 PAD 101PD Adapter*
Adult Radiotransparent, large posterior pad	A/B	P-214-M1		P-214-T1	P-214-Z1
Adult Radiotransparent pads and leadwires, large posterior pad, sterile	A/B	P-224-T1 PAD 1111 Adapter*	P-224-T1 PAD 109 Adapter*	P-224-T1	P-224-Z1
Pediatric Radiotranslucent pads, sterile	C	P-122-T1 PAD 1111 Adapter*	P-122-T1 PAD 109 Adapter*	P-122-T1	P-122-T1 PAD 101PD Adapter*
Pediatric Radiotranslucent pads	C	P-112-T1 PAD 1111 Adapter*	P-112-T1 PAD 109 Adapter*	P-112-T1	P-112-T1 PAD 101PD Adapter*
Pediatric Radiotransparent pads and leadwires	C	P-212-M1			P-212-Z1

***Adapters are reusable**

- 1 – M1 Medtronic models are comparable to Medtronic® Physio-Control QUIK-COMBO™ electrodes.
- 2 – P2 Philips models are comparable to Philips® gray flat-style double barrel electrodes.
- 3 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)
- 4 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

For Pediatric Use: ≥10kg (22lbs)
AED or AED Mode: ≥25kg (55lbs) ≥8 years old

Pad Dimensions:
A=6"x4.25" B=6"x9" C=4.25"x2.876"



TZ Medical
Sparked by your ideas

SAFETY & EFFICIENCY

TZ Medical's line of Safety and Efficiency products are designed to simplify and enhance staff and patient safety with new technology that improves ergonomics, comfort, and communication to provide an optimal user experience.

Pro11 Wireless Headset System • Anti-Fatigue Mats



PRO11 HEADSET SYSTEM

The TZ Medical Pro11 Headset System enables teams to talk hands-free in the Cath Lab and EP Lab and adjacent rooms - no raised voices or straining to hear, just clear, calm communication between every member of the team, reducing stress and creating a better working environment. The system has full connectivity with external broadcast and audiovisual streaming facilities, enabling the physician to participate real-time in remote broadcasts and seminars. Using high-quality DECT frequencies, multiple systems can be used in the same premises without interference.



Item Number	Description	UOM
-------------	-------------	-----

TZ-Q-P11BS	Pro11 Base Station	1 Each
TZ-Q-P11CH	Pro11 Charger	1 Each
TZ-Q-P11HS	Pro11 Headset (battery and headband included)	1 Each

Item Number	Description	UOM
-------------	-------------	-----

TZ-Q-P11BAT	Battery for Headset	1 Each
TZ-Q-P11HB	Headband for Headset	1 Each
TZ-Q-P11NB	Neckband for Headset	1 Each
TZ-Q-RH7HC	Hygiene Covers	Box of 100
TZ-Q-RH7EC	Ear Cushions	Box of 10
TZ-Q-P10SHELF	Battery Charger Shelf	1 Each




ANTI-FATIGUE MATS



TZ Medical
Sparked by your ideas


Features

- Unmatched ergonomic comfort & support
- Constructed with eco-friendly Bio-Foam®
- Durable, stain-resistant surface with brushed texture design
- Certified by the National Floor Safety Institute for superior high-traction bottom surface
- Edges beveled for added safety



Reduce injuries while improving productivity and safety

Available in 6 standard sizes as well as custom sizing



Beveled Edge
beveled at 18° to reduce tripping potential and makes rolling things on and off easy and safe

Matte Textured
top surface provides optimum friction



Durable
commercial grade heavy-duty polyurethane

Anti-microbial
properties for added protection

3/4" thick



High-Traction Bottom
certified by the National Floor Safety Institute

Integral Skin
stain and puncture resistant surface prevents mat from fluid absorption

Premium Comfort
proprietary Cellulon® core technology will not lose bounce over time

Eco-Friendly
constructed with Bio-Foam® a renewable, plant-based resource

Item Number	Description	UOM
NLEP-20X32	Durable 3/4" Thick Bio-Foam® Size 20X32, Black	1 Each
NLEP-20X48	Durable 3/4" Thick Bio-Foam® Size 20X48, Black	1 Each
NLEP-20X72	Durable 3/4" Thick Bio-Foam® Size 20X72, Black	1 Each
NLEP-24X36	Durable 3/4" Thick Bio-Foam® Size 24X36, Black	1 Each
NLEP-36X60	Durable 3/4" Thick Bio-Foam® Size 36X60, Black	1 Each
NLEP-36X114	Durable 3/4" Thick Bio-Foam® Size 36X114, Black	1 Each





TZ Medical
Sparked by your ideas

ELECTROPHYSIOLOGY


Working hand in hand with our cardiology specific products, TZ Medical has developed a variety of EP products that cover the spectrum of clinical EP applications.

EZ Holder • Booker Box • 14 Pole Catheter

EZ HOLDER



TZ Medical
Sparked by your ideas



The TZ Medical EZ Holder is a versatile device that acts as an anchor point to firmly hold catheters and lead wires in place during procedures. This innovative product will help reduce procedure time as well as improve the efficiency of EP cases and implants.

Holds your catheters, so you don't have to!

Features

- Helps prevent catheter movement during cardioversion, eliminating the need to re-position catheters after therapy is delivered
- Holds the catheter in place for HIS Bundle or any evoked potentials
- Frees up Physician's hands during ablation procedures
- Will help reduce chances of lead dislodgement during implant procedures
- Will hold a standard catheter or implant lead that is 5-9 French



Item Number	Description	UOM
CT-0010	Catheter and lead wire anchor	Box of 10



V1-MPC BOOKER BOX



The V1-MPC Patient Companion is a multi-port connector with switching capabilities for connection to defibrillators, EP study equipment and RF electrosurgical generators. The box provides a stable focal location for multiple wires, cables and connections.

The Booker Box System Provides:

- Multi-port Connector
- Patient Defibrillation
- Cardioversion
- Pacing... and more!



CALL US TODAY TO SET UP AN EVALUATION



Item Number	Description	UOM
-------------	-------------	-----


14 POLE CATHETER



TZ Medical
Sparked by your ideas

The dual-site 7F catheter provides the electrophysiologist with multiple options for the treatment of atrial and ventricular tachyarrhythmias. It allows the physician to position the distal end of the catheter into the coronary sinus (CS) while allowing the proximal electrodes to map the high right atrium (HRA).

Subclavian access eliminates the need for an additional femoral stick. The 14 platinum electrode spacing helps facilitate assessment of a line of conduction block.



Features

- Mapping, pacing and stimulation of atrial anatomy requiring large surface area electrodes
- Routine diagnostic EP procedures
- Pre-existing or induced atrial flutter and atrial fibrillation
- Clinical settings where arrhythmias are encountered

Diagnostic Features

Great handling reduces procedure time
Atraumatic introduction and removal
Maintains integrity for patient safety
High fidelity diagnostic electrograms

Unique Features

Fourteen large surface area electrodes
Dual-site configuration, 3 (70) 3 mm spacing
Mapping of both the CS and HRA
Quick connect cable

Item Number	Description	Configuration	mm	mm	Spacing	Type
-------------	-------------	---------------	----	----	---------	------

CF-D7-AAA-NSDA	Catheter-Map/Record/Pace	14-pole	4	4	3(70)3	Cournand Curve
----------------	--------------------------	---------	---	---	--------	----------------



TZ Medical
Sparked by your ideas

RADIATION PROTECTION

TZ Medical's line of radiation protection products provide excellent defense and peace of mind by protecting hospital staff from the harmful effects of secondary and scatter radiation.

Radiation Protection Eyewear ● Comfort Shield
Flipper Radiation Pad ● SR Headband



Barrier Lite

These polycarbonate wrap around frames weigh 60 grams, and are available with side shields. They have a lead equivalence of 0.75 mm Pb. Available in Crimson, Black, and Silver.



Atomic

These wrap-around glasses weigh just 60 grams, and have a lead equivalence of 0.75 mm Pb. The rubberized, anti-slip frames are antimicrobial, and feature vents to prevent fogging. Available in Black and Gunmetal Grey.



Astro II

The frames feature adjustable frame angle and temple length. Available with side shields for additional protection. The frames weigh 90 grams, and have a lead equivalence of 0.75 mm Pb. Available in Blue and Pink.



Oakley Crosslink

These sleek frames are constructed of light weight, durable O-Matter™ material, with Unobtainium™ nose pads and ear socks for enhanced grip. The frames weigh just 58 grams, and have a lead equivalence of 0.75 mm Pb. Available in Satin Black and Grey Smoke.



Nike Brazen

This lightweight and stylish design has a universally fitting frame for all face sizes, with a ventilated nose bridge for additional support. The frames weigh 76 grams, with a lead equivalence of 0.75 mm Pb. Available in Black and Tortoise.



Nike Rabid

This lightweight frames are designed to fit medium to large faces, with secure wrap temple arms for stability and comfort. The frames weigh 64 grams, with a lead equivalence of 0.75 mm Pb. Available in Black and Grey-Anthracite.



For a full list of Radiation Protection Eyewear and RX options, visit www.tzmedical.com



COMFORT SHIELD


The Comfort Shield is a leaded shield accessory that supplements the Comfort Cradle arm board. The shield measures 23" x 30", offering a large area of radiation protection to medical staff during procedures involving fluoroscopy.

The Comfort Shield is designed to fit easily over the TZ Medical Comfort Cradle arm board and attenuates nearly 98% of secondary and scatter radiation while remaining unobtrusive to the medical staff. It is easy to sanitize, and comes with a hanger for easy storage.



Features

- Attenuates 97.9% of secondary and scatter radiation at 80 kVp and has 0.5mm Lead Equivalence
- Low-profile design of Comfort Cradle maintained when shield is in use
- Easy to slip on and off the Comfort Cradle arm board
- Made of durable, flexible material that is easy to sanitize
- Comes with hanger for easy storage



Item Number	Description	UOM
-------------	-------------	-----

CZ-105

Comfort Shield


1 Each

FLIPPER



TZ Medical
Sparked by your ideas

The Flipper is a versatile radiation shield designed to keep physicians safe without getting in the way of the procedure. The Flipper's 19"x16" versatile design allows for use on either the right or left side for femoral groin access or subclavian access. Convenient adhesive strips on either side of the shield allow it to be secured in place to avoid slipping.



Lead free, and attenuating up to 75% of radiation at 90 kVp, the Flipper is ideal for Interventional Radiology, Cardiac Cath Labs, Electrophysiology and more

Features

- Attenuates up to 75% of radiation at 90 kVp*
- Reduces scatter radiation exposure to operators and staff
- Can be used for femoral groin access or subclavian access
- Adhesive on both sides to prevent slipping
- Lead equivalence 0.125mm
- Perfect for Interventional Radiology, Cardiac Cath, Electrophysiology & more



Item Number	Description	UOM
-------------	-------------	-----

CZ-007

Flipper Radiation Shield

Box of 15



SR HEADBAND

When it comes to radiation, protecting your brain is a top priority. The TZ Medical SR Headband is a cost effective way to block scatter radiation from reaching one of the most sensitive areas of the body during diagnostic and interventional procedures.

The SR Headband is disposable, lead free, easily worn, adjustable to any head size, and attenuates 87% of radiation at 80 kVp.* It's a must-have for physicians, techs, and anyone involved in imaging procedures.



Features

- Attenuates 87% of scatter radiation at 80 kVp*
- Comfortable and lightweight
- Lead free and disposable
- Fully adjustable, can be worn over any style of cap
- Cost effective way to increase safety in the lab



Item Number	Description	UOM
-------------	-------------	-----



PROCEDURE LIGHTING

TZ Medical offers a variety of lighting solutions to enhance visualization and improve procedure efficiency. Our line of high-quality examination lights provides multiple options for mobile, ceiling, and wall-mount systems, as well as single-use adhesive lights for poorly lit surgical cavities.

Wall Mounted Lights • Mobile Lights
Ceiling Mounted Lights • TruBrite



StarTrol LED Medical Lights produce a bright white light, (94 CRI and 4700K). They are ideal for true color rendition, and are cool and comfortable to work under, with virtually no radiant heat being generated by the LED's. The lights are energy-efficient, using less than one-third the energy of a halogen light, and require minimal maintenance. StarTrol LED Medical Lights are durable, precision machined in the USA and feature a 10-year warranty. TZ Medical offers a wide variety of light configurations, including ceiling-mounted, wall-mounted, and mobile rolling units, which are sure to fit into any workspace.

Features

- Energy Efficient & Eco-friendly
- Bright, natural white light
- 10-year limited warranty
- Cool to work under
- 40,000 hours of LED life
- Infinite variable positioning
- Multi-level dimming



Item Number	Description	UOM
ST-4X4-GSM	4 Light Spring Arm Mobile	Each
ST-8X4-GSM	8 Light Spring Arm Mobile	Each
ST-4X4-GSW	4 Light Spring Arm Wall	Each
ST-8X4-GSW	8 Light Spring Arm Wall	Each
ST-4X4-GS1	4 Light Spring Arm Single	Each
ST-8X4-GS1	8 Light Spring Arm Single	Each
ST-4X4-GS2	4 Light Spring Arm Dual	Each
ST-8X4-GS2	8 Light Spring Arm Dual	Each

Ceiling Mounted Units



TZ Medical
Sparked by your ideas



Dual Arm Ceiling Mounted Units



Wall Mounted Units



Mobile Units





Adhesive bottom for easy placement on various surfaces




Simple pull tab activates the light



Features

- Sterile, compact & disposable design that weighs less than 4 grams
- Provides continuous light for over 2 hours
- Adjustable LED light, providing optimal visibility for any procedure
- Self-adhesive for effortless mounting on multiple surfaces
- Ideal for pacemaker pockets with Army-Navy retractors



LED light can be positioned straight or bent for optimal illumination

Item Number	Description	UOM
-------------	-------------	-----

ARL-1010

Adhesive Retractor Light

Box of 10



TZ Medical
Sparked by your ideas

CARDIAC MONITORS

The TZ Medical line of cardiac monitoring devices combine wireless communication with easy-to-use monitors that help record and send important cardiac data in real time, helping physicians diagnose and treat their patients quickly and accurately.

Trident™ Pro • Trident™ Nano



TRIDENT PRO

Trident™ Pro cardiac monitors combine decades of research and innovation to deliver a solution designed to streamline remote ECG recording from start to finish. Built-in telecom and rechargeable lithium battery deliver the simplicity of a single-piece device with options for remote or local configuration and data transmission.



FEWER COMPONENTS. SMARTER MONITORING.

Trident™ Pro monitors are equipped with a large LED touch screen and on-board patient diary that provides a user-friendly experience for manually recorded patient events. Three channels of ECG data are acquired from surface electrodes adhered to the body, and 8G of internal memory with advanced detection algorithm offer a versatile solution to capture, store, and transmit every beat during Holter, Extended Holter, Event, and Mobile Cardiac Telemetry studies for up to 30 days.



Records up to 30 days of continuous EKG data



Large LCD color touchscreen with an on-board digital diary



Lithium battery may be charged while in use



Compact and discrete device with water resistant case



MONITORS	ITEM #	DESCRIPTION
Trident™ 40	H3R-0039	Digital Holter Monitor for Holter only studies
Trident™ 40L	H3R-0130	Digital 4-in-1 monitor for Holter, Extended Holter, MCT, and Event studies, LTE wireless coverage

TRIDENT NANO




TZ Medical
Sparked by your ideas

The Trident™ Nano is an accurate, reliable, water-resistant, and fully rechargeable monitor designed to capture clear and comprehensive ECG tracings during Holter and Extended Holter studies. 3 channels are recorded from either a lanyard cable or long-wear patch electrode, and its slim, unobtrusive design optimizes patient comfort and compliance.

With multiple wear options, 8G of internal storage, and a low-profile design, the Trident™ Nano offers a unique combination of minimal body burden with consistently clear ECG tracings, and the superior diagnostic yield of a 3-channel recording for study durations up to 21 days.

**For study durations lasting up to 21 days, multiple patch electrodes are needed.*



Records 3 channels of ECG from a single patch electrode



Fast and simple data retrieval via USB connection



Fully reusable and rechargeable device



3-color LED light indicator for device notifications



Water-resistant for showering

MONITORS	ITEM #	DESCRIPTION
Trident™ 30	HPR-0014	Digital Cardiac Monitor for Holter and Extended Holter studies, used with patch electrode and lanyard cable.



TZ Medical
Sparked by your ideas

CLEARBEAT MONITORING SERVICE

The Trident™ line of cardiac monitors are designed from the ground up to simplify and streamline the hookup process. Connect the device to USB, enter patient and study information into the ClearBeat web portal, then attach the device to the patient. It's that simple. When the device is returned, simply reconnect the USB to backup stored data and restore the device to default settings for future enrollments.

INNOVATION MEETS SIMPLICITY

CLEARBEAT SERVICE



TZ Medical
Sparked by your ideas



Fast & Reliable Arrhythmia Detection

Our certified cardiac technicians provide accurate and reliable analysis of Holter, Extended Holter, Event, and Mobile Cardiac Telemetry data 24/7, rapid notification of emergent events, and industry leading report turnaround times.



Customer Support

Our friendly customer service representatives are available to serve you and your patients 24/7. No lengthy teleprompt, no unnecessary hold-times. Just live patient and provider support when you need it, every time.



- Simplified enrollments
- One-step device configuration & reset
- Industry leading report turnaround times
- Review & e-sign reports from any web-connected device or facility EMR system
- Seamless transition from Holter to Event or MCT

**Improving cardiac telemetry.
Improving patients' lives.**



TZ Medical 
ClearBeat Monitoring Service

ITEM NUMBER INDEX

RADIAL ACCESS

ARC Radial Cuff

Box of 5
Item# RCB-0047 (Small)
Item# RCB-0048 (Medium)
Item# RCB-0049 (Large)

Cobra Board

Item# CBR-0004

StandTall™

Box of 5
Item# RD-ST0010L (10 cm)
Item# RD-ST0015L (15 cm)
Item# RD-ST0025L (25 cm)

Total Access Radial Board

Dimensions: 44.4" x 22.9"
Item# CZ-400

Secure Radial Board

Dimensions: 44.4" x 22.9"
Item# CZ-400-TVA

Disposable Straps
Dimensions: Two 2" x 24" straps,
One 2" x 18" strap
Item# TVA-MC-2

Radial Runway

Reusable Base
Item# RRU-0008
Disposable top pads, straps, and brace
Box of 10
Item# RRU-0009

FEMORAL ACCESS

Neptune Hemostasis Pad

Box of 10
Item# 8870-01
Neptune Plus
Box of 10
Item# 8870-03
SilverStat
Box of 10
Item# 8870-07

EZ Hold

Item# 8260-20

Pressure Discs

Box of 50
Item# 8880-01
Item# 8880-02
Item# 8880-03
Item# 8880-04
Item# 8880-06

Neptune Pressure Discs

Box of 10
Item# 8880-07
Item# 8880-08

PATIENT POSITIONING

Patient Immobilizer

Body Wrap
Item# CZ-106
GE Style Table
Item# CZ-107
Wrist Wraps
Box of 10
Item# CZ-008
XXL Extender
Item# CZ-012

PATIENT POSITIONING CONT.

Pannus Retention System

PRS-02
Box of 10
Item# CZ-PRS-02 (2 pads and 1 strap)
PRS-04
Box of 10
Item# CZ-PRS-04 (4 pads and 2 straps,
sterile)

Comfort Halo
Dimensions: 23" x 14.5"
Item# CZ-009

Sidewinder Board
Item# CZ-009-HH

Comfort Cradle
Dimensions: 21"x 13"x 9"
Item# CZ-100 (set of 2)
Item# CZ-001 (1 each)

Polypropylene Comfort Cradle
Dimensions: 21"x 13"x 9"
Item# CZ-100-PP (set of 2)
Item# CZ-001-PP (1 each)

Comfort Cloud
Box of 20
Item# CZ-002

SoftGels

Head Gels

Adult Gel Donut
8"OD x 3"ID x 1.75"H
Item# TZ-GELDA

Extra Deep Donut
8"OD x 3"ID x 3"H
Item# TZ-GELDA-3

Pediatric Donut
5.5"OD x 2.25"ID x 1.25"H
Item# TZ-GELPD

Neo Natal Donut
3.07"OD x 1.68"ID x .81"H
Item# TZ-GELND

Adult Horseshoe Donut
8"OD x 3"ID x 1.75"H
Item# TZ-GELHRS

Pediatric Horseshoe Donut
5.5"OD x 2.25"ID x 1.25"H
Item# TZ-GELPHR

Supine Headrest
11"OD x 9.5"ID x 2.625"H
Item# TZ-GELSHP

Head Pad
11"OD x 9.5"ID x 3.5"H
Item# TZ-GELHDP

Arm Gels & Ankle & Foot Gels
Flat Armboard Gel
24" x 4.5" x 0.625"
Item# TZ-GELAF

Gel Contoured Armboard
20" x 5.75" x 2.25"
Item# TZ-GELAC

Gel Ulnar Protector, Large
20.13" x 13.75" x .25"
Item# TZ-GELUPL

Gel Ulnar Protector, Small
16.25" x 6.5" x 0.5"
Item# TZ-GELUPS

PATIENT POSITIONING CONT.

SoftGels Cont.

Arm Gels & Ankle & Foot Gels Cont.

Gel Calf Cup
Item# TZ-GELCC

Gel Heel Cup
Item# TZ-GELHC

Heel Protector Pad
Item# TZ-GELHP

Ankle Stirrup Pad
Item# TZ-GELAS

Body Positioners

Gel Patient Positioner 10"
10" x 4.5" x 2" (Sandbag Style)
Item# TZ-GELPP10

Gel Patient Positioner 14"
14" x 4.5" x 2" (Sandbag Style)
Item# TZ-GELPP14

Gel Patient Positioner 18"
18" x 4.5" x 2" (Sandbag Style)
Item# TZ-GELPP18

Flat Bottom Gel Chest Roll 12"
12" x 4" x 4"
Item# TZ-GELCR4-12

Flat Bottom Gel Chest Roll 16"
16" x 4" x 4"
Item# TZ-GELCR4-16

Flat Bottom Gel Chest Roll 20"
20" x 4" x 4"
Item# TZ-GELCR4-20

Flat Bottom Gel Chest Roll 12"
12" x 6" x 6"
Item# TZ-GELCR6-12

Flat Bottom Gel Chest Roll 16"
16" x 6" x 6"
Item# TZ-GELCR6-16

Flat Bottom Gel Chest Roll 20"
20" x 6" x 6"
Item# TZ-GELCR6-20

Table Gels

Full Length Table Gel Overlay
20" x 72" x 0.5"
Item# TZ-GELFT

Gel Foot Table Overlay
20" x 20" x 0.5"
Item# TZ-GELFT3-FT

Gel Head Table Overlay
20" x 10" x 0.5"
Item# TZ-GELFT3-HD

Gel Torso Table Overlay
20" x 42" x 0.5"
Item# TZ-GELFT3-T

3 Piece Table Gel Overlay
Includes: Gel Foot Table Overlay, Gel Head
Table Overlay, Gel Torso Table Overlay
Item# TZ-GELFT3

INTERVENTIONAL

AccuFlate
Box of 5
Item# IDD-0001

M-CATH Microcatheter
Item# MCATH

INTERVENTIONAL CONT.

Control Mechanical Thrombectomy System

Item# ASP-1001

Item# ASP-1005

Item# ASP-1009

Item# ASP-1013

Item# ASP-1021

Control Mechanical Aspirator

Item# ASP-1029

Shoot 'N See - available in 3 quantities

Item# RCM-0005-10 (Box of 10)

Item# RCM-0005-30 (Box of 30)

Item# RCM-0005-50 (Box of 50)

ELECTRODES

Adult Electrodes

Radiotranslucent Electrode, leads out

Box of 10 (Size A)

Item# P-111-M1

Item# P-111-P2

Item# P-111-Z1

Item# P-111-T1

Radiotransparent Electrode

Box of 10 (Size A)

Item# P-211-M1

Item# P-211-P2

Item# P-211-Z1

Item# P-211-T1

Radiotransparent Electrode, leads out

Box of 10 (Size A)

Item# P-311-M1

Item# P-311-Z1

Radiotransparent Electrode, large back pad

Box of 10 (Size A & B)

Item# P-214-M1

Item# P-214-Z1

Item# P-214-T1

Pediatric Electrodes

Radiotranslucent Electrode

Box of 5 (Size C)

Item# P-112-T1

Radiotransparent Electrode

Box of 5 (Size C)

Item# P-212-M1

Item# P-212-Z1

Sterile Electrodes

Sterile Radiotranslucent Electrode

Box of 10

Item# P-121-T1

Sterile Radiotransparent Electrode w/ large back pad

Box of 10

Item# P-224-T1

Item# P-224-Z1

Sterile Radiotranslucent Pediatric Electrode

Box of 10

Item# P-122-T1

SAFETY & EFFICIENCY

Pro11 Headset

Pro11 Base Station

Item# TZ-Q-P11BS

Pro11 Charger

Item# TZ-Q-P11CH

Pro11 Headset (w/Battery & Headband)

Item# TZ-Q-P11HS

SAFETY & EFFICIENCY

Pro11 Headset Cont.

Battery for Headset

Item# TZ-Q-P11BAT

Headband for Headset

Item# TZ-Q-P11HB

Hygiene Covers

Box of 100

Item# TZ-Q-RH7HC

Ear Cushions

Box of 10

Item# TZ-Q-RH7EC

Battery Charger Shelf

Item# TZ-Q-P10SHELF

Anti-Fatigue Mats

Durable ¾" Thick, Black Bio-foam®,

Size 20x32

Item# NLEP-20X32

Durable ¾" Thick, Black Bio-foam®,

Size 20x48

Item# NLEP-20X48

Durable ¾" Thick, Black Bio-foam®,

Size 20x72

Item# NLEP-20X72

Durable ¾" Thick, Black Bio-foam®,

Size 24x36

Item# NLEP-24X36

Durable ¾" Thick, Black Bio-foam®,

Size 36x60

Item# NLEP-36X60

Durable ¾" Thick, Black Bio-foam®,

Size 36x114

Item# NLEP-36X114

ELECTROPHYSIOLOGY

EZ Holder

Box of 10

Item# CT-0010

Booker Box

Item# V1-MPC

Catheter

Item# CF-D7-AAA-NSDA

RADIATION PROTECTION

Eyewear

Astro II

Item# AST-II

Atomic

Item# ATOM

Circuit XL

Item# CIRCL

Barrier Lite

Item# BLT

MX-30

Item# MX

Costa Brine

Item# BRINE

Nike Brazen

Item# NBRA

Nike Rabid

Item# NRAB

RADIATION PROTECTION CONT.

Eyewear Cont.

Oakley Crosslink

Item# XLINK

Oakley Jupiter Squared

Item# JUPITER

Oakley Straightlink

Item# SLINK

Rudy Project Contender

Item# CONTENDER

Comfort Shield

Item# CZ-105

Flipper

Box of 15

Item# CZ-007

SR Headband

Box of 10

Item# CZ-013

PROCEDURE LIGHTING

StarTrol

StarTrol 4X4-GSW

Wall Examination Light

Item# ST-4X4-GSW

StarTrol 8X4-GSW

Wall Examination Light

Item# ST-8X4-GSW

StarTrol 4X4-GS1

Ceiling Mounted Light

Item# ST-4X4-GS1

StarTrol 8X4-GS1

Ceiling Mounted Light

Item# ST-8X4-GS1

StarTrol 4X4-GS2

Ceiling Mounted Light

Item# ST-4X4-GS2

StarTrol 8X4-GS2

Ceiling Mounted Light

Item# ST-8X4-GS2

StarTrol 4X4-GSM

Mobile Examination Light

Item# ST-4X4-GSM

StarTrol 8X4-GSM

Mobile Examination Light

Item# ST-8X4-GSM

TruBrite

Box of 10

Item# ARL-1010

CARDIAC MONITORS

Trident™ Pro

Trident™ 40

Holter Monitor

Item# H3R-0039

Trident™ 40L

4-1 Device - Holter, Extended Holter, MCT, and Event Monitor

Item# H3R-0130

Trident™ 30

Holter and Extended Holter Monitor

Item# HPR-0014

For additional information or questions

about ordering call: 800.944.0187

or visit us at

www.tzmedical.com

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

15396 R14.4



TZ Medical
Sparked by your ideas