

Trident™ Pro

Trident™ Pro cardiac monitors combine decades of research and innovation to deliver a solution designed to streamline remote ECG recording from start to finish. Both our dedicated Holter-only and 4-in-1 recorders with built-in LTE cellular modem deliver the simplicity and ease-of-use of a single piece monitoring solution.



MONITORS	ITEM #	DESCRIPTION
Trident™ 40	H3R-0039	Digital Holter monitor for Holter & Extended Holter studies only
Trident™ 40L	H3R-0130	4-in-1 Device - Holter, Extended Holter, MCT, and Event Monitor

Trident™ Pro monitors are equipped with a large LED touch screen and on-board patient diary that provides a user-friendly experience for manually recorded patient events. Three channels of ECG data are acquired from surface electrodes adhered to the body, and 8G of internal memory with advanced detection algorithm offer a versatile solution to capture, store, and transmit every beat for up to 30 days.


3-Electrode Cable & Charger



3-Electrode Cable with Charging Port & Barrel Jack Charger (Charge-While-Wear)



9-Bay Charger



USB Data Cable

Belt Clip Holster

Fewer Components. Smarter Monitoring.



Records up to 30 days of continuous EKG data



Large LCD color touchscreen with an on-board digital diary



Rechargeable Lithium Battery



Compact and discrete device with water resistant case

Accessories