



TZ Medical
Sparked by your ideas

URPRO HANGER

The URPro Hanger provides a stable platform for mounting URPro Adapter Panels beneath imaging table controls, helping protect staff from scatter radiation while maintaining normal table operation. Available in two sizes, it offers a lightweight, modular approach to lower-body shielding and allows damaged panels to be replaced individually rather than replacing an entire shield assembly.



- Attenuates 98% of secondary and scatter radiation at 80 kVp
- Lead equivalent of 0.50 mm
- Attaches directly to table-side rails and positions shielding beneath table controls
- Allows normal operation and access to table controls during procedures
- Supports multiple modular URPro panels for customizable protection
- Lightweight design ideal for Cath Labs, IR suites, OBLs, and ASCs

Shield Staff Without Blocking Table Controls



Designed for Everyday Protection

Position shielding where staff need it most while maintaining full access to table controls.

Item #	Description	UOM
BDR-1091	URPro Hanger - 4 Panel, 30" L	1 Each
BDR-1092	URPro Hanger - 5 Panel, 30" L	1 Each
BDR-1093	URPro Hanger - 4 Panel, 35" L	1 Each
BDR-1094	URPro Hanger - 5 Panel, 35" L	1 Each