




TZ Medical
Sparked by your ideas


URPRO OVERHEAD



The TZ Medical URPro Overhead is a universal shield that directly replaces your current MAVIG™ overhead shielding with high-quality, long-lasting radiation protection. The URPro Overhead shield is well-structured with fabric binding around the edges, and flexible overlapping strips with a lead equivalency of 0.50 mm.

Features

- Attenuates 98% of secondary and scatter radiation at 80 kVp
- Lead equivalent of 0.50 mm
- Available in 40 cm and 50 cm widths
- Directly compatible with current MAVIG™ overhead shielding systems
- Simple installation process



Item Number	Description	UOM
BDR-1030	40cm Overhead URPro Panel Assembly	EACH
BDR-1035	50cm Overhead URPro Panel Assembly	EACH

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com
24248 R0.3