




TZ Medical
Sparked by your ideas

URPRO UPPER SHIELD PANELS



Enhance your radiation protection setup with TZ Medical's URPro Upper Shield Replacement Panels, engineered for full compatibility with the MAVIG™ UT69 framing system. Crafted from high-quality, long-lasting materials, each panel is built to withstand daily clinical use while maintaining exceptional protective performance.



These panels feature a durable fabric binding around every edge for added structural integrity and a clean, professional finish. With a lead equivalency of 0.50 mm, the URPro panels deliver reliable attenuation to help safeguard healthcare professionals during interventional procedures.

Features

- Attenuates 98% of secondary and scatter radiation at 80 kVp
- Lead equivalent of 0.50 mm
- Directly compatible with current MAVIG™ UT69 shielding systems
- High-quality and long-lasting
- Available in small (17cm) and large (25cm) sizes



Item Number	Description	UOM
BDR-1053	URPro Upper Shield Panels, Small	Set of 4
BDR-1052	URPro Upper Shield Panels, Large	Set of 4

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

00000 R0.1